



CONEFOR

CONEFOR is an electrosurgical electrode for the cervical conization, called cone biopsy.

CONEFOR – OBTAIN PRECISE CONE BIOPSIES EFFORTLESSLY AND SAFELY.

The freely rotating hub of the CONEFOR with three prongs will stabilize the tissue during the biopsy.

The diameter and height of the **CONEFOR** will determine the biopsy size.

CONEFOR Biopsy Dimensions

Height: 9, 12, 15, 18, 21 mm

Diameter: 12, 18, 24, 30 mm

Length of Instrument: 175 mm

Shaft diameter: 3/32"

Packed 5/box - weight of the box:

150 grams - dimension of the box:

100 x 50 x 290mm

Conization Electrode

CONEFOR			
Models			
No.	Size	Dia.	Height
231812	Small	18mm	12mm
231815	Small – extended	18mm	15mm
231818	Small - deep	18mm	18mm
232412	Medium	24mm	12mm
232415	Medium – extended	24mm	15mm
232418	Medium - deep	24mm	18mm
233009	Large	30mm	9mm
233012	Large – extended	30mm	12mm
233015	Large – deep	30mm	15mm
231221	Mini -- extended	12mm	21mm

New England Medical

Phone 845-778-4200

Distributed by: