




ELECTROSURGICAL ACCESSORIES

TECHNICAL INFORMATION

INDEX

Page 06:	<i>Power Supply Cable</i>
Page 09:	<i>Surgical pencil</i>
Page 12:	<i>Foot Switch</i>
Page 16:	<i>Neutral</i>
Page 23:	<i>Connection</i>
Page 29:	<i>Electrode</i>
Page 32:	<i>Adapters</i>
Page 34:	<i>Bipolar/Monopolar</i>
Page 41:	<i>Evacuation Kit</i>
Page 45:	<i>Units support</i>

CAPTIONS

Compatibility:

+ = optional

x = not compatible

= to verify

N = Quantity Supplied



Packaging:
quantity of articles contained



Autoclavable at 134°



Not reusable



Sterilized with ethylene

IP67

**The item is completely dust-proof
and can with stand a 30-minute
immersion in standing water**



Anti explosion



100 Sterilization cycles



30 Sterilization cycles





ACCESSORIES



POWER SUPPLY CABLE


Description Power supply cable SIE-IEC
Cable Lenght 5 m (196,85 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC



pcs: 1

00100.01




Description Power supply cable SIE-IEC
Cable Lenght 2 m (78,74 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC



pcs: 1

00100.03




Description Power supply cable IT-IEC
Cable Lenght 2 m (78,74 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC



pcs: 1

00100.00




Description Power supply cable USA-IEC
Cable Lenght 2 m (78,74 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC




pcs: 1

00100.04



00100.05




Description Power supply cable UK-IEC
Cable Lenght 2 m (78,74 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC



pcs: 1

00100.07




Description Power supply cable BR-IEC
Cable Lenght 2 m (78,74 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC



pcs: 1

00100.09




Description Power supply cable AU/CN-IEC
Cable Lenght 2 m (78,74 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC



pcs: 1

00100.10




Description Power supply cable JP-IEC
Cable Lenght 2 m (78,74 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC




pcs: 1



SURGICAL PENCIL

00201.02



Description	Handle for microsurgical needle autoclavable
Cable Lenght	2 m (78,74 in)
Size	10x1x1 cm (3,93x0,39x0,39 in)
Weight	0,05 kg (0,11lb)
Electrode Compatibility	ø0,8 mm (0,031 in)
Sterilization	Autoclavable
Sterilization cycles	100
Compatibility	500500.L



pcs: 1

00205.00



Description	Handle with switches
Cable Length	3 m (118,11 in)
Size	18,5x1,6x1,7 cm (7,2x0,62x0,66 in)
Weight	0,03kg (0,06 lb)
Electrode Compatibility	ø2,4 mm (0,94 in)
Sterilization	Autoclavable
Sterilization cycles	100



pcs: 1

00205.40



Description	Handle with switches
Cable Length	3 m (118,11 in)
Size	18,5x1,6x1,7 cm (7,2x0,62x0,66 in)
Weight	0,03kg (0,06 lb)
Electrode Compatibility	ø4 mm (1,15 in)
Sterilization	Autoclavable
Sterilization cycles	100



pcs: 1

Description Monopolar Handpiece
Cable Length 3 m (118,11 in)
Size 18,5x1,6x1,7 cm (7,2x0,62x0,66 in)
Weight 0,03kg (0,06 lb)
Electrode Compatibility ø2,4 mm (0,94 in)
Sterilization Autoclavable
Sterilization cycles 100



pcs: 1

50 D
1
80 D
+
120
+
160
+
200
+
TOUCH 200
+
FLASH 120
+
FLASH 160
+
FLASH 200
+
300 HP
+
400 HP
+
TOUCH HP
+
TOUCH HP (Argon)
x
SURTRON EVAC
x
SURTRON ABC

00206.00



Description Monopolar Handpiece
Cable Length 3m (118,11 in)
Size 18,5x1,6x1,7 cm (7,2x0,62x0,66 in)
Weight 0,03 kg (0,06 lb)
Electrode Compatibility ø4 mm (1,15 in)
Sterilization Autoclavable
Sterilization cycles 100



pcs: 1

50 D
+
80 D
+
120
+
160
+
200
+
TOUCH 200
+
FLASH 120
+
FLASH 160
+
FLASH 200
+
300 HP
+
400 HP
+
TOUCH HP
+
TOUCH HP (Argon)
x
SURTRON EVAC
x
SURTRON ABC

00206.40 S




Description Handpiece for argon electrodes reusable
Cable Length 3,5 m (137,8 in)
Size 18,5x1,6x1,7 cm (7,2x0,62x0,66 in)
Weight 0,03kg (0,06lb)
Electrode Compatibility ø4 mm (1,15 in)
Sterilization Autoclavable
Sterilization cycles 100



pcs: 1

50 D
x
80 D
x
120
x
160
x
200
x
TOUCH 200
x
FLASH 120
x
FLASH 160
x
FLASH 200
x
300 HP
x
400 HP
x
TOUCH HP
x
TOUCH HP (Argon)
+
SURTRON EVAC
+
SURTRON ABC

00420.00



Description Disposable handle with finger switches (F4797)
Size 14x2x1 cm (5,5x0,78x0,39 in)
Cable Length 3m (118,11in)
Weight 0,07 kg (0,15 lb)
Electrode Compatibility ø2,4 mm (0,094 in)
Sterilization Sterile EO
Sterilization cycles 0



pcs: 1

50 D
x
80 D
+
120
+
160
+
200
+
TOUCH 200
+
FLASH 120
+
FLASH 160
+
FLASH 200
+
300 HP
+
400 HP
+
TOUCH HP
+
TOUCH HP (Argon)
x
SURTRON EVAC
x
SURTRON ABC

755VL





FOOT SWITCH

Description
Cable Length
Size
Weight

Water-proof foot switch
 2 m (78,75 in)
 11x3,5x9 cm (7,2x0,62x0,66 in)
 0,3 kg (0,66 lb)

50 D
 1
 80 D
 1
 80
 1
 120
 1
 160
 1
 200
 1
 TOUCH 200
 1
 FLASH 120
 1
 FLASH 160
 1
 FLASH 200
 1
 300 HP
 x
 400 HP
 x
 TOUCH HP
 x
 TOUCH HP (Argon)
 x
 SURTRON EVAC
 +
 SURTRON ABC
 x



pcs: 1

00304.00



WATER PROOF IP67

Description
Cable Length
Size
Weight


Single water-proof foot switch
 (3 pins)
 5 m (196,85 in)
 11x3,5x9 cm (7,2x0,62x0,66 in)
 0,3 kg (0,66 lb)

50 D
 x
 80 D
 x
 80
 x
 120
 x
 160
 x
 200
 x
 TOUCH 200
 x
 FLASH 120
 x
 FLASH 160
 x
 FLASH 200
 x
 # 300 HP
 x
 400 HP
 #
 TOUCH HP
 x
 TOUCH HP (Argon)
 x
 SURTRON EVAC
 x
 SURTRON ABC
 +



pcs: 1

00304.01



WATER PROOF IP67

Description
Cable Length
Size
Weight


Water-proof foot switch
 3 m (118,11 in)
 11x3,5x9 cm (7,2x0,62x0,66 in)
 0,3 kg (0,66 lb)

+ 50 D
 + 80 D
 + 80
 + 120
 + 160
 + 200
 + TOUCH 200
 + FLASH 120
 + FLASH 160
 + FLASH 200
 + 300 HP
 x
 400 HP
 x
 TOUCH HP
 x
 TOUCH HP (Argon)
 x
 SURTRON EVAC
 +
 SURTRON ABC
 x




pcs: 1

00304.04



WATER PROOF IP67

00304.PO



Description Protection for single foot switch
Cable Length NO
Size 12x10x14 cm (4,7x3,9x1,64 in)
Weight 0,65 kg (0,143 lb)

50 D +
80 D +
80 +
120 +
160 +
200 +
TOUCH 200 +
FLASH 120 +
FLASH 160 +
FLASH 200 +
300 HP x
400 HP x
TOUCH HP x
TOUCH HP (Argon) x
SURTRON EVAC #
SURTRON ABC #



pcs: 1

00305.03



Description Double water proof foot switch
Cable Length 5 m (196,85 in)
Size 35x6x18 cm (13,78x2,36x7 in)
Weight 1.3 kg (2,86lb)

50 D x
80 D x
80 +
120 +
160 +
200 +
TOUCH 200 +
FLASH 120 +
FLASH 160 +
FLASH 200 +
300 HP x
400 HP x
TOUCH HP x
TOUCH HP (Argon) x
SURTRON EVAC x
SURTRON ABC x



pcs: 1

00305.04



Description Waterproof double foot switch
Cable Length 5 m (196,85 in)
Size 35x6x18 cm (13,78x2,36x7 in)
Weight 1.3 kg (2,86lb)

50 D x
80 D x
80 x
120 x
160 x
200 x
TOUCH 200 x
FLASH 120 x
FLASH 160 x
FLASH 200 x
300 HP 2
400 HP 2
TOUCH HP x
TOUCH HP (Argon) x
SURTRON EVAC x
SURTRON ABC x



pcs: 1

WATER PROOF IP67



**NEU
TRAL**


Description	Rod neutral electrode with cable
Cable Length	1,8 m (70,86 in)
Size	12xø1,8 cm (4,72x0,7 in)
Weight	0,1 kg (0,22 lb)
Sterilization	NO
Sterilization cycles	NO

+ 50 D	x TOUCH 200
+ 80 D	x FLASH 120
+ 80	x FLASH 160
+ 120	x FLASH 200
+ 160	x 300 HP
x 200	x 400 HP
x TOUCH 200	x TOUCH HP
x TOUCH HP	x TOUCH HP (Argon)
x SURTRON EVAC	x SURTRON EVAC
x SURTRON ABC	x SURTRON ABC



pcs: 1

00400.00



MAX. 120W

Integrated connection


Description	Steel neutral electrode with cable
Cable Length	3 m (118,11 in)
Size	16x12 cm (6,3x4,7 in)
Weight	0,1 kg (0,22 lb)
Sterilization	NO
Sterilization cycles	NO

+ 50 D	x TOUCH 200
+ 80 D	x FLASH 120
+ 80	x FLASH 160
+ 120	x FLASH 200
+ 160	x 300 HP
+ 200	x 400 HP
x TOUCH 200	x TOUCH HP
x TOUCH HP	x TOUCH HP (Argon)
x SURTRON EVAC	x SURTRON EVAC
x SURTRON ABC	x SURTRON ABC



pcs: 1

00401.00



Integrated connection


Description	Steel neutral electrode with cable
Cable Length	3 m (118,11 in)
Size	24x16 cm (9,4x6,3 in)
Weight	0,1 kg (0,22 lb)
Sterilization	NO
Sterilization cycles	NO

+ 50 D	x TOUCH 200
+ 80 D	x FLASH 120
+ 80	x FLASH 160
+ 120	x FLASH 200
+ 160	x 300 HP
+ 200	x 400 HP
x TOUCH 200	x TOUCH HP
x TOUCH HP	x TOUCH HP (Argon)
x SURTRON EVAC	x SURTRON EVAC
x SURTRON ABC	x SURTRON ABC




pcs: 1

00401.01



Integrated connection

00401.02



Integrated connection


Description	Steel neutral electrode autoclavable
Size WxHxD	16x12x0,4 cm (6,3x4,7x0,15 in)
Cable Length	3 m (118,11 in)
Weight	0,1 kg (0,22 lb)
Sterilization	Autoclavable
Sterilization cycles	100

50 D	+	TOUCH 200
80 D	+	FLASH 120
80	+	FLASH 160
120	+	FLASH 200
160	+	300 HP
200	+	400 HP
	+	TOUCH HP
	+	TOUCH HP (Argon)
	+	SURTRON EVAC
	+	SURTRON ABC



pcs: 1

00401.03



Integrated connection

Description	Steel neutral electrode autoclavable
Size WxHxD	24x16x0,4 cm (9,4x6,3x0,15 in)
Cable Length	3 m (118,11 in)
Weight	0,1 kg (0,22 lb)
Sterilization	Autoclavable
Sterilization cycles	100

50 D	+	TOUCH 200
80 D	+	FLASH 120
80	+	FLASH 160
120	+	FLASH 200
160	+	300 HP
200	+	400 HP
	+	TOUCH HP
	+	TOUCH HP (Argon)
	+	SURTRON EVAC
	+	SURTRON ABC



pcs: 1

00401.10



Description	Steel neutral electrode FLEX
Size WxHxD	12x21x0,2 cm (4,7x8,2x0,078 in)
Cable Length	NO
Weight	0,01 Kg (0,02lb)
Sterilization	Autoclavable
Sterilization cycles	100
Compatibility	00404.08, 00404.09, 00404.10, 00404.11

50 D	+	TOUCH 200
80 D	+	FLASH 120
80	+	FLASH 160
120	+	FLASH 200
160	+	300 HP
200	+	400 HP
	+	TOUCH HP
	+	TOUCH HP (Argon)
	+	SURTRON EVAC
	+	SURTRON ABC




pcs: 1

Description	Steel neutral electrode FLEX with cable
Size WxHxD	12x21x0,2 cm (4,7x8,2x0,078 in)
Cable Length	3m (118,11 in)
Weight	0,01 Kg (0,02lb)
Sterilization	NO
Sterilization cycles	NO

+ 50 D	+ TOUCH 200
+ 80 D	+ FLASH 120
+ 80	+ FLASH 160
+ 120	+ FLASH 200
+ 160	+ 300 HP
+ 200	+ 400 HP
+ TOUCH 200	+ TOUCH HP
+ FLASH 120	+ TOUCH HP (Argon)
+ FLASH 160	x SURTRON EVAC
+ FLASH 200	x SURTRON ABC
+ 300 HP	
+ 400 HP	
+ TOUCH HP	
+ TOUCH HP (Argon)	
x SURTRON EVAC	
x SURTRON ABC	



00401.11



Integrated connection

Description	Steel neutral electrode FLEX with cable autoclavable
Size WxHxD	12x21x0,2 cm (4,7x8,2x0,078 in)
Cable Length	3m (118,11 in)
Weight	0,01 Kg (0,02lb)
Sterilization	Autoclavable
Sterilization cycles	100

+ 50 D	+ TOUCH 200
+ 80 D	+ FLASH 120
+ 80	+ FLASH 160
+ 120	+ FLASH 200
+ 160	+ 300 HP
+ 200	+ 400 HP
+ TOUCH 200	+ TOUCH HP
+ FLASH 120	+ TOUCH HP (Argon)
+ FLASH 160	x SURTRON EVAC
+ FLASH 200	x SURTRON ABC
+ 300 HP	
+ 400 HP	
+ TOUCH HP	
+ TOUCH HP (Argon)	
x SURTRON EVAC	
x SURTRON ABC	



00401.12



Integrated connection

134°C




Description	Steel neutral electrode FLEX S
Size WxHxD	12x21x0,2 cm (4,7x8,2x0,078 in)
Cable Length	NO
Weight	0,01 Kg (0,02 lb)
Sterilization	Autoclavable
Sterilization cycles	100
Compatibility	00404.08, 00404.09, 00404.10, 00404.11

+ 50 D	+ TOUCH 200
+ 80 D	+ FLASH 120
+ 80	+ FLASH 160
+ 120	+ FLASH 200
+ 160	+ 300 HP
+ 200	+ 400 HP
+ TOUCH 200	+ TOUCH HP
+ FLASH 120	+ TOUCH HP (Argon)
+ FLASH 160	x SURTRON EVAC
+ FLASH 200	x SURTRON ABC
+ 300 HP	
+ 400 HP	
+ TOUCH HP	
+ TOUCH HP (Argon)	
x SURTRON EVAC	
x SURTRON ABC	




00401.20



134°C



00401.21



Integrated connection

Description

Steel neutral electrode FLEX S with cable

Size WxHxD

12x21x0,2 cm (4,7x8,2x0,078 in)

Cable Length

3 m (118,11 in)

Weight

0,01 Kg (0,02 lb)

Sterilization

NO

Sterilization cycles


NO

#	50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
	+	+	+	+	+	+	+	+	+	+	+	+	+	+	x	x



pcs: 1

00401.22



Integrated connection



Description

Steel neutral electrode FLEX S with cable autoclavable

Size WxHxD

12x21x0,2 cm (4,7x8,2x0,078 in)

Cable Length

3 m (118,11 in)

Weight

0,01 kg (0,02 lb)

Sterilization

Autoclavable

Sterilization cycles


100

#	50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
	+	+	+	+	+	+	+	+	+	+	+	+	+	+	x	x



pcs: 1

00403.01



Description

Steel neutral electrode autoclavable

Size WxHxD

10x1,2x1,2 cm (3,934x0,47x0,47 in)

Cable Length

NO

Weight

0,1 kg (0,22 lb)

Sterilization

Autoclavable

Sterilization cycles

100

Compatibility

00404.08, 00404.09, 00404.10, 00404.11

1	50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x



pcs: 1

Description	Rod neutral electrode
Size WxHxD	14x2x2 cm (5,5x0,78x0,78 in)
Cable Lenght	NO
Weight	0,4 Kg (0,88 lb)
Sterilization	Autoclavable
Sterilization cycles	100
Compatibility	00404.02, 00404.04

+ 50 D	
x 80 D	
x 80	
x 120	
x 160	
x 200	
x TOUCH 200	
x FLASH 120	
x FLASH 160	
x FLASH 200	
x 300 HP	
x 400 HP	
x TOUCH HP	
x TOUCH HP (Argon)	
x SURTRON EVAC	
x SURTRON ABC	



pcs: 1

00403.02




Description	Steel neutral electrode
Size WxHxD	16x12x0,02 cm (6,3x4,7x0,008 in)
Cable Lenght	NO
Weight	0,1 kg (0,22 lb)
Sterilization	Autoclavable
Sterilization cycles	100
Compatibility	00404.08, 00404.09, 00404.10, 00404.11

+ 50 D	
1 80 D	
1 80	
1 120	
1 160	
1 200	
1 TOUCH 200	
1 FLASH 120	
1 FLASH 160	
1 FLASH 200	
x 300 HP	
x 400 HP	
x TOUCH HP	
x TOUCH HP (Argon)	
x SURTRON EVAC	
x SURTRON ABC	



pcs: 1

5365A




Description	Steel neutral electrode
Size WxHxD	24x16x0,02 cm (9,4x6,3x0,008 in)
Cable Lenght	NO
Weight	0,1 kg (0,22 lb)
Sterilization	Autoclavable
Sterilization cycles	100
Compatibility	00404.08, 00404.09, 00404.10, 00404.11

+ 50 D	
1 80 D	
1 80	
1 120	
1 160	
1 200	
1 TOUCH 200	
1 FLASH 120	
1 FLASH 160	
1 FLASH 200	
x 300 HP	
x 400 HP	
+ TOUCH HP	
+ TOUCH HP (Argon)	
x SURTRON EVAC	
x SURTRON ABC	




pcs: 1

6429A



F7920




Description	Disposable split neutral electrode (F7820)
Size WxHxD	20x10,5x0,2 cm (7,8x4,1x0,078 in)
Cable Lenght	NO
Weight	0,01 Kg (0,02 lb)
Sterilization	NO
Sterilization cycles	NO
Compatibility	00404.08, 00404.09, 00404.10, 00404.11

#	50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
	+	+	+	+	+	+	+	+	+	+	+	+	+	+	x	x



pcs: 1

F7930



Description	Conductive rubber split neutral electrode
Size WxHxD	20x10,5x0,4 cm (7,8x4,1x0,15 in)
Cable Lenght	NO
Weight	0,1 Kg (0,2lb)
Sterilization	NO
Sterilization cycles	NO
Compatibility	00404.07

#	50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
	+	+	+	+	+	+	+	+	+	+	+	+	+	+	x	x




pcs: 1




CON NECT ION


00402.00




00402.01



00402.02



00402.03



LD 24

Description	Monopolar Cable M4-F2.8
Weight	0,06 kg (0,13lb)
Cable Lenght	3 m (118,11 in)
Sterilization	Autoclavable
Sterilization cycles	100



pcs: 1

Description	Monopolar Cable M4-F4
Weight	0,06 kg (0,13 lb)
Cable Length	3 m (118,11 in)
Sterilization	Autoclavable
Sterilization cycles	100



pcs: 1

Description	Monopolar Cable M4-MP4
Weight	0,06 kg (0,13 lb)
Cable Length	3 m (118,11 in)
Sterilization	Autoclavable
Sterilization cycles	100



pcs: 1

Description	Monopolar Cable M4-EU
Weight	0,06 kg (0,13 lb)
Cable Length	3 m (118,11 in)
Sterilization	Autoclavable
Sterilization cycles	100




pcs: 1

Description Monopolar Cable M4-F2.M
Weight 0,06 kg (0,13 lb)
Cable Length 3 m (118,11 in)
Sterilization Autoclavable
Sterilization cycles 100

50 D
 # 80 D
 # 80
 # 120
 # 160
 # 200
 + TOUCH 200
 + FLASH 120
 + FLASH 160
 + FLASH 200
 + 300 HP
 + 400 HP
 + TOUCH HP
 + TOUCH HP (Argon)
 x SURTRON EVAC
 x SURTRON ABC



pcs: 1




Description Monopolar Cable M4-F2.M
Weight 0,06 kg (0,13 lb)
Cable Length 1,8 m (70,86 in)
Sterilization NO
Sterilization cycles NO

50 D
 1 x 80 D
 x 80
 x 120
 x 160
 x 200
 x TOUCH 200
 x FLASH 120
 x FLASH 160
 x FLASH 200
 x 300 HP
 x 400 HP
 x TOUCH HP
 x TOUCH HP (Argon)
 x SURTRON EVAC
 x SURTRON ABC



pcs: 1




Description Cable for rod neutral electrode autoclavable
Cable Length 3 m (118 in)
Unit Connection 4 mm (0,157 in)
Sterilization Autoclavable
Sterilization cycles 100

50 D
 + x 80 D
 x 80
 x 120
 x 160
 x 200
 x TOUCH 200
 x FLASH 120
 x FLASH 160
 x FLASH 200
 x 300 HP
 x 400 HP
 x TOUCH HP
 x TOUCH HP (Argon)
 x SURTRON EVAC
 x SURTRON ABC



pcs: 1




Description Cable for neutral plate F7930
Cable Length 5m (196,85 in)
Unit Connection 4 mm (0,157 in)
Accessory Connection 2,8 mm (0,110 in)
Compatibility F7930

50 D
 + 80 D
 + 80
 + 120
 + 160
 + 200
 + TOUCH 200
 + FLASH 120
 + FLASH 160
 + FLASH 200
 + 300 HP
 + 400 HP
 + TOUCH HP
 + TOUCH HP (Argon)
 x SURTRON EVAC
 x SURTRON ABC



pcs: 1




00402.04

00404.02

00404.04

00404.07

00404.08




Description	Cable for neutral electrode disposable type 5365/6429/FLEX
Weight	0,06 kg (0,13 lb)
Cable Lenght	3 m (118,11 in)
Sterilization	NO
Sterilization cycles	NO
Compatibility	0350, F790, 00401.10, 00401.20, 5365A, 6429A

+ 50 D	TOUCH 200	TOUCH HP (Argon)
+ 1 80 D	FLASH 120	SURTRON EVAC
+ 1 80	FLASH 160	SURTRON ABC
+ 1 120	FLASH 200	
+ 1 160	300 HP	
+ 1 200	400 HP	
+ 1 TOUCH 200	TOUCH HP	
+ 1 FLASH 120	TOUCH HP (Argon)	
+ 1 FLASH 160	SURTRON EVAC	
+ 1 FLASH 200	SURTRON ABC	
+ 1 300 HP		
+ 1 400 HP		
+ 1 TOUCH HP		
+ 1 TOUCH HP (Argon)		
+ 1 SURTRON EVAC		
+ 1 SURTRON ABC		



pcs: 1

00404.09




Description	Cable for connection neutral disposable 5365-6429/FLEX FLEX S AUTOCLAVE
Weight	0,06 kg (0,13 lb)
Cable Lenght	3 m (118,11 in)
Sterilization	Autoclavable
Sterilization cycles	100
Compatibility	0350, F790, 00401.10, 00401.20, 5365A, 6429A

+ 50 D	TOUCH 200	TOUCH HP (Argon)
+ 1 80 D	FLASH 120	SURTRON EVAC
+ 1 80	FLASH 160	SURTRON ABC
+ 1 120	FLASH 200	
+ 1 160	300 HP	
+ 1 200	400 HP	
+ 1 TOUCH 200	TOUCH HP	
+ 1 FLASH 120	TOUCH HP (Argon)	
+ 1 FLASH 160	SURTRON EVAC	
+ 1 FLASH 200	SURTRON ABC	
+ 1 300 HP		
+ 1 400 HP		
+ 1 TOUCH HP		
+ 1 TOUCH HP (Argon)		
+ 1 SURTRON EVAC		
+ 1 SURTRON ABC		



pcs: 1

00404.10




Description	Cable for neutral plate US type
Weight	0,06 kg (0,13 lb)
Cable Lenght	5 m (169,85 in)
Sterilization	NO
Sterilization cycles	NO
Compatibility	0350, F790, 00401.10, 00401.20, 5365A, 6429A

50 D	TOUCH 200	TOUCH HP (Argon)
# 80 D	FLASH 120	SURTRON EVAC
# 80	FLASH 160	SURTRON ABC
# 120	FLASH 200	
# 160	300 HP	
# 200	400 HP	
# TOUCH 200	TOUCH HP	
# FLASH 120	TOUCH HP (Argon)	
# FLASH 160	SURTRON EVAC	
# FLASH 200	SURTRON ABC	
# 300 HP		
# 400 HP		
# TOUCH HP		
# TOUCH HP (Argon)		
# SURTRON EVAC		
# SURTRON ABC		



pcs: 1

00404.11



Description	Cable for neutral plate US type AUTOCLAVE
Weight	0,06 kg (0,13 lb)
Cable Lenght	5 m (169,85 in)
Sterilization	Autoclavable
Sterilization cycles	100
Compatibility	0350, F790, 00401.10, 00401.20, 5365A, 6429A

50 D	TOUCH 200	TOUCH HP (Argon)
# 80 D	FLASH 120	SURTRON EVAC
# 80	FLASH 160	SURTRON ABC
# 120	FLASH 200	
# 160	300 HP	
# 200	400 HP	
# TOUCH 200	TOUCH HP	
# FLASH 120	TOUCH HP (Argon)	
# FLASH 160	SURTRON EVAC	
# FLASH 200	SURTRON ABC	
# 300 HP		
# 400 HP		
# TOUCH HP		
# TOUCH HP (Argon)		
# SURTRON EVAC		
# SURTRON ABC		



pcs: 1


Description Bipolar cable EUR
Weight 0,06 kg (0,13 lb)
Cable Length 3 m (118,11 in)
Sterilization NO
Sterilization cycles NO
Compatibility 310 Bipolar Forceps series

50 D
x
80 D
x
80
+
120
+
160
+
200
+
TOUCH 200
+
FLASH 120
+
FLASH 160
+
FLASH 200
+
300 HP
+
400 HP
+
TOUCH HP
+
TOUCH HP (Argon)
x
SURTRON EVAC
x
SURTRON ABC



pcs: 1

00411.00



Description Bipolar cable TWIN
Weight 0,06 kg (0,13 lb)
Cable Length 3 m (118,11 in)
Sterilization NO
Sterilization cycles NO
Compatibility 310 Bipolar Forceps series

50 D
x
80 D
x
80

120

160

200

TOUCH 200
x
FLASH 120

FLASH 160


FLASH 200

300 HP
x
400 HP
x
TOUCH HP
x
TOUCH HP (Argon)
x
SURTRON EVAC
x
SURTRON ABC



pcs: 1

00412.00



Description Artery Sealer Bipolar cable
Weight 0,06 kg (0,13 lb)
Cable Length 3 m (118,11 in)
Sterilization NO
Sterilization cycles NO
Compatibility 110-750NS, 110-755NS, 110-760NS,

50 D
x
80 D
x
80

120

160

200

TOUCH 200
+
FLASH 120


FLASH 160

FLASH 200
+
300 HP
+
400 HP
+
TOUCH HP
+
TOUCH HP (Argon)
x
SURTRON EVAC
x
SURTRON ABC



pcs: 1

00413.00



Description Bipolar cable US 3m
Weight 0,60 kg (0,13 lb)
Cable Length 3 m (118,11 in)
Sterilization NO
Sterilization cycles NO
Compatibility 110-700

50 D
x
80 D
x
80

120

160

200

TOUCH 200

FLASH 120

FLASH 160

FLASH 200

300 HP

400 HP


TOUCH HP

TOUCH HP (Argon)
x
SURTRON EVAC
x
SURTRON ABC




pcs: 1

00414.00



00415.00




Description Bipolar cable ENDO (MP3 - F3)
Weight 0,06 kg (0,13 lb)
Cable Length 3 m (118,11 in)
Sterilization Autoclavable
Sterilization cycles 100

50 D x
 80 D x
 80 x
 120 x
 160 x
 200 +
 TOUCH 200 +
 FLASH 120 x
 FLASH 160 +
 FLASH 200 +
 300 HP +
 400 HP +
 TOUCH HP +
 TOUCH HP (Argon) +
 SURTRON EVAC x
 SURTRON ABC x



pcs: 1

00416.00



Description Bipolar cable ENDO (MP3 - F3)
Weight 0,06 kg (0,13 lb)
Cable Length 3 m (118,11 in)
Sterilization Autoclavable
Sterilization cycles 100

50 D x
 80 D x
 80 x
 120 x
 160 x
 200 +
 TOUCH 200 +
 FLASH 120 x
 FLASH 160 +
 FLASH 200 +
 300 HP +
 400 HP +
 TOUCH HP +
 TOUCH HP (Argon) +
 SURTRON EVAC x
 SURTRON ABC x



pcs: 1

00417.00




Description Bipolar cable ENDO3 (MP2-F4)
Weight 0,06 kg (0,13 lb)
Cable Length 3 m (118,11 in)
Sterilization Autoclavable
Sterilization cycles 100

50 D x
 80 D x
 80 x
 120 x
 160 x
 200 +
 TOUCH 200 +
 FLASH 120 x
 FLASH 160 +
 FLASH 200 +
 300 HP +
 400 HP +
 TOUCH HP +
 TOUCH HP (Argon) +
 SURTRON EVAC x
 SURTRON ABC x



pcs: 1

00418.00



Description Bipolar cable ENDO3 (MP2-F4)
Weight 0,06 kg (0,13 lb)
Cable Length 3 m (118,11 in)
Sterilization Autoclavable
Sterilization cycles 100

50 D x
 80 D x
 80 x
 120 x
 160 x
 200 +
 TOUCH 200 +
 FLASH 120 x
 FLASH 160 +
 FLASH 200 +
 300 HP +
 400 HP +
 TOUCH HP +
 TOUCH HP (Argon) +
 SURTRON EVAC x
 SURTRON ABC x




pcs: 1




ELEC TRODE


00500.00



500500.L11



152 Short



Description Kit of assorted electrode length 5 cm (10pcs)

SINGLE ITEM

- 500500.L1: Straight thin wire electrode (5pcs)
- 500500.L3: Loop electrode $\varnothing 4$ mm (5pcs)
- 500500.L4: Loop electrode $\varnothing 6$ mm (5pcs)
- 500500.L7: Drop electrode length (5pcs)
- 500500.L8: Noose electrode length (5pcs)
- 500500.L9: Straight ball electrode $\varnothing 3$ mm (5pcs)
- 500500.L2: Bent thin wire electrode (5pcs)
- 500500.L5: Drop electrode length (5pcs)
- 500500.L6: Bent thick wire electrode (5pcs)
- 500500.L10: Bent5 ball electrode $\varnothing 3$ mm (5pcs)

Compatibility 00205.00, 00206.00, 755VL



kit 00500.00: assorted pcs 10
for each items from L1 to L10: pcs 5

50 D	80 D	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
1	1	1	1	1	1	1	1	1	x	+	+	+	x	x

Description Needles for micro-surgery length 1cm
Weight 0,002 kg (0,004 lb)
Size 12x0,05x0,8 mm (0,47x0,002x0,03 in)
Sterilization Autoclavable
Sterilisation cycles NO

50 D	80 D	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
+	+	+	+	+	+	1	1	1	x	x	#	#	#	#



pcs: 10

Description

152-110 Blade Electrode 7 cm

50 D	80 D	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP
+	+	+	+	+	+	+	+	+	+	+	#

- 152-112 Blade Curved Electrode 7 cm
- 152-120 Needle Curved Electrode 7 cm
- 152-122 Needle Electrode 13 cm
- 152-130 Ball Electrode $\varnothing 2$ mm 6 cm
- 152-132 Ball Curved Electrode $\varnothing 2$ mm 5 cm
- 152-140 Ball Electrode $\varnothing 3$ mm 6cm
- 152-142 Ball Curved Electrode $\varnothing 3$ mm 6 cm

50 D	80 D	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP
+	+	+	+	+	+	+	+	+	+	+	#

152-152 Ball Curved Electrode $\varnothing 4$ mm 6 cm

50 D	80 D	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP
x	+	+	+	+	+	+	+	+	+	+	#

- 152-160 Ball Electrode $\varnothing 5$ mm 6 cm
- 152-162 Ball Curved Electrode $\varnothing 5$ mm 6cm

50 D	80 D	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP
x	x	+	+	+	+	+	+	+	+	+	#

Compatibility 00205.00, 00206.00, 755VL

Description Kit of assorted electrode length 10 cm (10pcs)

SINGLE ITEM

500500.L/L1: Straight thin wire electrode length (5pcs)

500500.L/L3: Loop electrode $\varnothing 4$ mm length (5pcs)

500500.L/L4: Loop electrode $\varnothing 6$ mm (5pcs)

500500.L/L7: Drop electrode length (5pcs)

500500.L/L8: Noose electrode length (5pcs)

500500.L/L9: Straight ball electrode $\varnothing 3$ mm (5pcs)

500500.L/L2: Bent thin wire electrode (5pcs)

500500.L/L5: Bent hook electrode (5pcs)

500500.L/L6: Bent thick wire electrode (5pcs)

500500.L/L10: Straight ball electrode $\varnothing 3$ mm (5pcs)

Compatibility 00205.00, 00206.00, 755VL


+ 50 D
+ 80 D
+ 80
+ 120
+ 160
+ 200
+ TOUCH 200
+ FLASH 120
+ FLASH 160
+ FLASH 200
+ 300 HP
+ 400 HP
+ TOUCH HP
+ TOUCH HP (Argon)
+ SURTRON EVAC
+ SURTRON ABC



pcs:

kit 00500.00_L/L: assorted pcs 10
for each items from L-L1 to L-L10 pcs 5

00500.00/L



Description

152-115 Blade Electrode 16 cm

152-125 Needle Electrode 13 cm

152-145 Ball Electrode $\varnothing 3$ mm 13 cm

152-165 Ball Electrode $\varnothing 5$ mm 14 cm

152-175-10 Loop Electrode 10x10 l.14 cm

152-190-13 Loop Electrode 20x13 l.14 cm

152-190-20 Loop Electrode 20x20 l.15 cm

152-195 Conization Electrode 13 cm


Compatibility 00205.00, 00206.00, 755VL

+ 50 D
+ 80 D
+ 80
+ 120
+ 160
+ 200
+ TOUCH 200
+ FLASH 120
+ FLASH 160
+ FLASH 200
+ 300 HP
+ 400 HP
+ TOUCH HP
+ TOUCH HP (Argon)
+ SURTRON EVAC
+ SURTRON ABC



pcs: 5

152 Long





ADAP TERS

Description Bipolar adapter
Weight 0,03 kg (0,06 lb)
Size 4,5x1,5x7,5 cm (1,8x0,6x3 in)
Sterilization NO
Sterilization cycles NO
Compatibility 00414.00, 00418.00

50 D	TOUCH 200	TOUCH HP	TOUCH HP (Argon)
x x 80 D	x x FLASH 120	x x 400 HP	x x SURTRON EVAC
x x 80	x x FLASH 160	x x 300 HP	x x SURTRON ABC
x x 120	x x FLASH 200		
x x 160			
x x 200			



pcs: 1

00498.00



Description 2001 bipolar adapter selector (until 2009)
Weight 0,03 kg (0,06 lb)
Size 3x3x8,5 cm (1,2x1,2x3 in)
Sterilization NO
Sterilization cycles NO
Compatibility 00412.00

50 D	TOUCH 200	TOUCH HP	TOUCH HP Argon
x x 80 D	x x FLASH 120	x x 400 HP	x x SURTRON EVAC
x x 80	x x FLASH 160	x x 300 HP	x x SURTRON ABC
x x 120	x x FLASH 200		
x x 160			
x x 200			



pcs: 1

00498.04



Description Adapter for neutral electrode 6,3mm/Valley
Weight 0,015kg (0,03lb)
Size 2,2x2,2x7,2 cm (0,86x0,86x2,8 in)
Sterilization NO
Sterilization cycles NO
Compatibility 00404.10, 00404.11

50 D	TOUCH 200	TOUCH HP	TOUCH HP (Argon)
+ 80 D	+ FLASH 120	+ 400 HP	+ SURTRON EVAC
+ 80	+ FLASH 160	+ 300 HP	+ SURTRON ABC
+ 120	+ FLASH 200		
+ 160			
+ 200			



pcs: 1

00498.06



Description Bipolar adapter EUR/2xM2,5
Weight 0,015 kg (0,03 lb)
Size 2,2x2,2x7,2 cm (0,86x0,86x2,8 in)
Sterilization NO
Sterilization cycles NO

50 D	TOUCH 200	TOUCH HP	TOUCH HP (Argon)
x x 80 D	x x FLASH 120	x x 400 HP	x x SURTRON EVAC
x x 80	x x FLASH 160	x x 300 HP	x x SURTRON ABC
x x 120	x x FLASH 200		
x x 160			
x x 200			



pcs: 1

00498.08




Description Bipolar adapter EUR/3xM4
Weight 0,015 kg (0,03 lb)
Size 2,2x2,2x7,2 cm (0,86x0,86x2,8 in)
Sterilization NO
Sterilization cycles NO

50 D	TOUCH 200	TOUCH HP	TOUCH HP (Argon)
# 80 D	# FLASH 120	# 400 HP	# SURTRON EVAC
# 80	# FLASH 160	# 300 HP	# SURTRON ABC
# 120	# FLASH 200		
# 160			
# 200			



pcs: 1

00498.10





BI
PO
LAR

MO
NO
PO
LAR

Description US Bi-Turbinate Electrode
Size 18,5x2x2 cm (7,3x0,78x78 in)
Weight 0,07 kg (0,154 in)
Tip Thickness 4,5 mm (1,77 in)
Sterilization Autoclavable
Sterilization cycles 30
Compatibility 00414.00

50 D x
 80 D x
 80 x
 120 x
 160 +
 200 x
 TOUCH 200 x
 FLASH 120 x
 FLASH 160 x
 FLASH 200 x
 300 HP x
 400 HP x
 TOUCH HP #
 TOUCH HP (Argon) #
 SURTRON EVAC #
 SURTRON ABC #



pcs: 1



175 mm



Description Bipolar Forceps Artery Sealer
Size 27 cm (10,62 in)
Weight 0,1 kg (0,22 lb)
Sterilization Autoclavable
Sterilization cycles 30
Tip thickness 3 mm
Compatibility 00413.00

50 D x
 80 D x
 80 x
 120 x
 160 +
 200 x
 TOUCH 200 x
 FLASH 120 x
 FLASH 160 x
 FLASH 200 x
 300 HP x
 400 HP x
 TOUCH HP #
 TOUCH HP (Argon) #
 SURTRON EVAC #
 SURTRON ABC #



pcs: 1



27 cm



Tip Shape



Description Bipolar Forceps Artery Sealer
Size 25,5 cm (10 in)
Size 27 cm (10,62 in)
Sterilization Autoclavable
Sterilization cycles 30
Tip thickness 3mm
Compatibility 00413.00

50 D x
 80 D x
 80 x
 120 x
 160 +
 200 x
 TOUCH 200 x
 FLASH 120 x
 FLASH 160 x
 FLASH 200 x
 300 HP x
 400 HP x
 TOUCH HP #
 TOUCH HP (Argon) #
 SURTRON EVAC #
 SURTRON ABC #



pcs: 1



25,5 cm



Tip Shape



Description Bipolar Forceps Artery Sealer
Size 17 cm (6,7in)
Size 27 cm (10,62 in)
Sterilization Autoclavable
Sterilization cycles 30
Tip thickness 2 mm
Compatibility 00413.00

50 D x
 80 D x
 80 x
 120 x
 160 +
 200 x
 TOUCH 200 x
 FLASH 120 x
 FLASH 160 x
 FLASH 200 x
 300 HP x
 400 HP x
 TOUCH HP #
 TOUCH HP (Argon) #
 SURTRON EVAC #
 SURTRON ABC #



pcs: 1



17 cm



Tip Shape



110-700

110-750NS

110-755NS

110-760NS

310-110-05



11,5 cm



Tip Shape



Description	Bipolar Forceps 11,5cm
Weight	0,03 kg (0,6lb)
Length	11,5 cm (4,52 in)
Tip thickness	0,5 mm (0,019in)
Sterilization	Autoclavable
Sterilization cycles	30
Compatibility	00411.00, 00412.00

x 50 D	
x 80 D	
x 80	
x 120	
x 160	
x 200	
+ TOUCH 200	
+ FLASH 120	
+ FLASH 160	
+ FLASH 200	
+ 300 HP	
+ 400 HP	
# TOUCH HP	
# TOUCH HP (Argon)	
# SURTRON EVAC	
# SURTRON ABC	



pcs: 1

310-112-05



11,5 mm



Tip Shape



Description	Curved Bipolar Forceps
Weight	0,03 kg (0,6lb)
Length	11,5 cm (4,52 in)
Tip thickness	0,5 mm (0,019in)
Sterilization	Autoclavable
Sterilization cycles	30
Compatibility	00411.00, 00412.00

x 50 D																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																										
--------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--



pcs: 1

310-140-10



20 cm



Tip Shape



Description	Bipolar Forceps
Weight	0,03 kg (0,6lb)
Length	20 cm (7,87lb)
Tip thickness	1 mm (0,039in)
Sterilization	Autoclavable
Sterilization cycles	30
Compatibility	00411.00, 00412.00

x 50 D	
x 80 D	
+ 80	
+ 120	
+ 160	
+ 200	
+ TOUCH 200	
+ FLASH 120	
+ FLASH 160	
+ FLASH 200	
+ 300 HP	
+ 400 HP	
# TOUCH HP	
# TOUCH HP (Argon)	
# SURTRON EVAC	
# SURTRON ABC	



pcs: 1

310-140-20



20 cm



Tip Shape



Description	Bipolar Forcep
Weight	0,03 kg (0,6lb)
Length	20 cm (7,87 lb)
Tip thickness	2 mm (0,078 in)
Sterilization	Autoclavable
Sterilization cycles	30
Compatibility	00411.00, 00412.00

x 50 D	
x 80 D	
+ 80	
+ 120	
+ 160	
+ 200	
+ TOUCH 200	
+ FLASH 120	
+ FLASH 160	
+ FLASH 200	
+ 300 HP	
+ 400 HP	
# TOUCH HP	
# TOUCH HP (Argon)	
# SURTRON EVAC	
# SURTRON ABC	



pcs: 1

Description Curved Bipolar Forceps
Weight 0,03Kg (0,6lb)
Length 20 cm (7,87 in)
Tip thickness 1 mm (0,039 in)
Sterilization Autoclavable
Sterilization cycles 30
Compatibility 00411.00, 00412.00

x 50 D
x 80 D
x 80
x 120
x 160
x 200
+ TOUCH 200
x FLASH 120
x FLASH 160
x FLASH 200
+ 300 HP
+ 400 HP
TOUCH HP
TOUCH HP (Argon)
SURTRON EVAC
SURTRON ABC



pcs: 1

310-142-10



20 cm



Tip Shape



Description Curved Bipolar Forceps
Weight 0,03 Kg (0,6 lb)
Length 20 cm (7,87 in)
Tip thickness 2 mm (0,078 in)
Sterilization Autoclavable
Sterilization cycles 30
Compatibility 00411.00, 00412.00

x 50 D
x 80 D
+ 80
+ 120
+ 160
+ 200
+ TOUCH 200
+ FLASH 120
+ FLASH 160
+ FLASH 200
+ 300 HP
+ 400 HP
TOUCH HP
TOUCH HP (Argon)
SURTRON EVAC
SURTRON ABC



pcs: 1

310-142-20



20 cm



Tip Shape



Description Bipolar Forceps Artery Sealer
Weight 0,03 kg (0,6 lb)
Length 25 cm (10 in)
Tip thickness 3 mm (0,11 in)
Sterilization Autoclavable
Sterilization cycles 30
Compatibility 00411.00, 00412.00

x 50 D
x 80 D
+ 80
+ 120
+ 160
+ 200
+ TOUCH 200
+ FLASH 120
+ FLASH 160
+ FLASH 200
+ 300 HP
+ 400 HP
TOUCH HP
TOUCH HP (Argon)
SURTRON EVAC
SURTRON ABC



pcs: 1

310-180-10



25 cm



Tip Shape



Description Angled Bipolar Forceps 20cm
Weight 0,03 kg (0,6 lb)
Length 17 cm (0,07 in)
Tip thickness 2 mm (0,07 in)
Sterilization Autoclavable
Sterilization cycles 30
Compatibility 00411.00, 00412.00

x 50 D
x 80 D
+ 80
+ 120
+ 160
+ 200
+ TOUCH 200
+ FLASH 120
+ FLASH 160
+ FLASH 200
+ 300 HP
+ 400 HP
TOUCH HP
TOUCH HP (Argon)
SURTRON EVAC
SURTRON ABC



pcs: 1

310-180-20



17 cm



Tip Shape



310-182-10



20 cm



310-185-10



20 cm



310-510



20 cm



310-550



20 cm



LD 38

Description	Angled Curved Bipolar Forceps
Weight	0,03 kg (0,6 lb)
Length	20 cm (7,87 lb)
Tip thickness	1 mm (0,039 in)
Sterilization	Autoclavable
Sterilisation cycles	30
Compatibility	00411.00, 00412.00



Description	Angled Curved Bipolar Forceps
Weight	0,03 kg (0,6lb)
Length	20 cm (7,87 lb)
Tip thickness	1 mm (0,039in)
Sterilization	Autoclavable
Sterilisation cycles	30
Compatibility	00411.00, 00412.00



pcs: 1

Description	Bipolar Electrode
Weight	0,03 kg (0,6 lb)
Length	20 cm (7,87 lb)
Tip thickness	1 mm (0,039 in)
Sterilization	Autoclavable
Sterilisation cycles	30
Compatibility	00411.00, 00412.00



pcs: 1

Description	Bipolar Electrode
Weight	0,03 kg (0,06lb)
Length	20 cm (7,87lb)
Tip thickness	0,5 mm (0,02in)
Sterilization	Autoclavable
Sterilisation cycles	30
Compatibility	00411.00, 00412.00



pcs: 1

Description Bipolar Forcep 20 cm - curved 2
Weight 0,03 Kg (0,6 lb)
Length 18 cm (7 lb)
Tip thickness 0,5 mm (0,02 in)
Sterilization Autoclavable
Sterilisation cycles 30
Compatibility 00411.00, 00412.00

x 50 D
x 80 D
+ 80
+ 120
+ 160
+ 200
+ TOUCH 200
+ FLASH 120
+ FLASH 160
+ FLASH 200
+ 300 HP
+ 400 HP
TOUCH HP
TOUCH HP (Argon)
SURTRON EVAC
SURTRON ABC



pcs: 1

310-590



18 cm

Tip Shape



Description Bipolar Forcep Clamp Scissors
Weight 0,05 kg (0,1 lb)
Length 15 cm (5,9 lb)
Tip thickness 0,5 mm (0,02 in)
Sterilization Autoclavable
Sterilisation cycles 30
Compatibility 00418.00

x 50 D
x 80 D
+ 80
+ 120
+ 160
+ 200
+ TOUCH 200
+ FLASH 120
+ FLASH 160
+ FLASH 200
+ 300 HP
+ 400 HP
TOUCH HP
TOUCH HP (Argon)
SURTRON EVAC
SURTRON ABC



pcs: 1

410-100-15



15 cm

Tip Shape



Description Bipolar Forcep Clamp Scissors
Weight 0,05 kg (0,1 lb)
Length 19 cm (7,5 lb)
Tip thickness 0,5 mm (0,02 in)
Sterilization Autoclavable
Sterilisation cycles 30
Compatibility 00418.00

x 50 D
x 80 D
+ 80
+ 120
+ 160
+ 200
+ TOUCH 200
+ FLASH 120
+ FLASH 160
+ FLASH 200
+ 300 HP
+ 400 HP
TOUCH HP
TOUCH HP (Argon)
SURTRON EVAC
SURTRON ABC



pcs: 1

410-100-19



19 cm

Tip Shape



Description Bipolar Forcep Scissors
Weight 0,05 kg (0,11 lb)
Length 18 cm (7,0 lb)
Tip thickness 2 mm (0,078 in)
Sterilization Autoclavable
Sterilisation cycles 30
Compatibility 00418.00

x 50 D
x 80 D
+ 80
+ 120
+ 160
+ 200
+ TOUCH 200
+ FLASH 120
+ FLASH 160
+ FLASH 200
+ 300 HP
+ 400 HP
TOUCH HP
TOUCH HP (Argon)
SURTRON EVAC
SURTRON ABC



pcs: 1

410-200-18




18 cm

Tip Shape



410-200-21



21 cm



Description	Bipolar Forcep Scissors
Weight	0,05 kg (0,11 lb)
Length	21 cm (7,0 lb)
Tip thickness	2 mm (0,078 in)
Sterilization	Autoclavable
Sterilisation cycles	30
Compatibility	00418.00

50 D	x	80 D	x	80	x	120	+	160	+	200	+	TOUCH 200	+	FLASH 120	+	FLASH 160	+	FLASH 200	+	300 HP	+	400 HP	#	TOUCH HP	#	TOUCH HP (Argon)	#	SURTRON EVAC	#	SURTRON ABC	#
------	---	------	---	----	---	-----	---	-----	---	-----	---	-----------	---	-----------	---	-----------	---	-----------	---	--------	---	--------	---	----------	---	------------------	---	--------------	---	-------------	---



pcs: 1

410-200-23



23 cm




Description	Bipolar Forcep Scissors
Weight	0,05 kg (0,11 lb)
Length	23 cm (7,0 lb)
Tip thickness	2 mm (0,078 in)
Sterilization	Autoclavable
Sterilisation cycles	30
Compatibility	00418.00

50 D	x	80 D	x	80	x	120	+	160	+	200	+	TOUCH 200	+	FLASH 120	+	FLASH 160	+	FLASH 200	+	300 HP	+	400 HP	#	TOUCH HP	#	TOUCH HP (Argon)	#	SURTRON EVAC	#	SURTRON ABC	#
------	---	------	---	----	---	-----	---	-----	---	-----	---	-----------	---	-----------	---	-----------	---	-----------	---	--------	---	--------	---	----------	---	------------------	---	--------------	---	-------------	---




pcs: 1

330-160



18 cm



Description	Monopolar Forceps 21cm
Weight	0,05Kg (0,11lb)
Length	18 cm (7,0 lb)
Cutter length	45 mm (1,7 in)
Sterilization	Autoclavable
Sterilisation cycles	30
Compatibility	0040202

50 D	x	80 D	x	80	x	120	+	160	+	200	+	TOUCH 200	+	FLASH 120	+	FLASH 160	+	FLASH 200	+	300 HP	+	400 HP	#	TOUCH HP	#	TOUCH HP (Argon)	#	SURTRON EVAC	#	SURTRON ABC	#
------	---	------	---	----	---	-----	---	-----	---	-----	---	-----------	---	-----------	---	-----------	---	-----------	---	--------	---	--------	---	----------	---	------------------	---	--------------	---	-------------	---




pcs: 1

330-134-20



21 cm



Description	Monopolar Surgical Scissors
Weight	0,05 kg (0,11 lb)
Length	21 cm (8,2 lb)
Tip thickness	2 mm (0,078 in)
Sterilization	Autoclavable
Sterilisation cycles	30
Compatibility	0040202

50 D	x	80 D	x	80	x	120	+	160	+	200	+	TOUCH 200	+	FLASH 120	+	FLASH 160	+	FLASH 200	+	300 HP	+	400 HP	#	TOUCH HP	#	TOUCH HP (Argon)	#	SURTRON EVAC	#	SURTRON ABC	#
------	---	------	---	----	---	-----	---	-----	---	-----	---	-----------	---	-----------	---	-----------	---	-----------	---	--------	---	--------	---	----------	---	------------------	---	--------------	---	-------------	---



pcs: 1



EVA CUA TION KIT

00900.01/12



STERILE EO 2

Description	Handle evacuation kit for surgical smoke evacuato
Pipe Diameter	ø22 cm (3,93 in)
Pipe Lenght	18,5 cm (7,2 in)
Weight	0,05 kg (0,11 lb)
Sterilization	NO
Sterilisation cycles	NO

50 D	x	80 D	x	80	x	120	x	160	x	200	x	TOUCH 200	x	FLASH 120	x	FLASH 160	x	FLASH 200	x	300 HP	x	400 HP	x	TOUCH HP	x	TOUCH HP (Argon)	x	SURTRON EVAC	+	SURTRON ABC	x
------	---	------	---	----	---	-----	---	-----	---	-----	---	-----------	---	-----------	---	-----------	---	-----------	---	--------	---	--------	---	----------	---	------------------	---	--------------	---	-------------	---



pcs: 12

00900.02/06



2

Description	Evacuation kit for surgical smoke evacuator
Pipe Diameter	ø22 cm (0,86 in)
Pipe Lenght	240 cm (94,48 in)
Weight	0,05 kg (0,11 lb)
Sterilization	NO
Sterilisation cycles	NO
Compatibility	00900.01/12, 00900.05, 00900.07, 00900.08/12,

50 D	x	80 D	x	80	x	120	x	160	x	200	x	TOUCH 200	x	FLASH 120	x	FLASH 160	x	FLASH 200	x	300 HP	x	400 HP	x	TOUCH HP	x	TOUCH HP (Argon)	x	SURTRON EVAC	+	SURTRON ABC	x
------	---	------	---	----	---	-----	---	-----	---	-----	---	-----------	---	-----------	---	-----------	---	-----------	---	--------	---	--------	---	----------	---	------------------	---	--------------	---	-------------	---



pcs: 6

00900.04



134°C 100 CYCLES

Description	Autoclavable evacuation tube D.22
Pipe Diameter	ø7 mm (0,27 in)
Pipe Lenght	2,40 m (94,5 in)
Weight	0,05 kg (0,11 lb)
Sterilization	Autoclavable
Sterilisation cycles	100
Compatibility	00900.01/12, 00900.05, 00900.06/06, 00900.07, 00900.08/12

50 D	x	80 D	x	80	x	120	x	160	x	200	x	TOUCH 200	x	FLASH 120	x	FLASH 160	x	FLASH 200	x	300 HP	x	400 HP	x	TOUCH HP	x	TOUCH HP (Argon)	x	SURTRON EVAC	+	SURTRON ABC	x
------	---	------	---	----	---	-----	---	-----	---	-----	---	-----------	---	-----------	---	-----------	---	-----------	---	--------	---	--------	---	----------	---	------------------	---	--------------	---	-------------	---



pcs: 1

Description Autoclavable evacuation tube D.22
Pipe Diameter ø7 mm (0,27 in)
Pipe Length 2,40 m (94,5 in)
Weight 0,05 kg (0,11lb)
Sterilization Autoclavable
Sterilisation cycles 100
Compatibility 00900.02, 00900.04

50 D	x	TOUCH 200	x
80 D	x	FLASH 120	x
80	x	FLASH 160	x
120	x	FLASH 200	x
160	x	300 HP	x
200	x	400 HP	x
	x	TOUCH HP	x
	x	TOUCH HP (Argon)	x
	1	SURTRON EVAC	
	x	SURTRON ABC	x



pcs: 1

00900.05



Description Bacterial/Viral Filter D.22
Pipe Diameter ø22 mm (0,867 in)
Connection Diameter ø15 mm (0,59 in)
Sterilization NO
Sterilisation cycles NO
Compatibility 00900.04, 00900.02/06, 00900.FU01/10

50 D	x	TOUCH 200	x
80 D	x	FLASH 120	x
80	x	FLASH 160	x
120	x	FLASH 200	x
160	x	300 HP	x
200	x	400 HP	x
	x	TOUCH HP	x
	x	TOUCH HP (Argon)	x
	+	SURTRON EVAC	
	x	SURTRON ABC	x



pcs: 1

00900.06/06



Description Autoclavable handle evacuation kit
Pipe Diameter ø22 cm (0,86 in)
Pipe Length 150 cm (59 in)
Weight 0,05 kg (0,11 lb)
Sterilization Autoclavable
Sterilisation cycles 30
Compatibility 755VL, 00205.00/40, 00206.00/40, 00900.02/06, 00900.04

50 D	x	TOUCH 200	x
80 D	x	FLASH 120	x
80	x	FLASH 160	x
120	x	FLASH 200	x
160	x	300 HP	x
200	x	400 HP	x
	x	TOUCH HP	x
	x	TOUCH HP (Argon)	x
	+	SURTRON EVAC	
	x	SURTRON ABC	x



pcs: 1

00900.07



00900.08/12



STERILE EO



Description	Evacuation kit for TROCAR Luer Lock
Pipe Diameter	ø22 mm (0,86 in)
Pipe Lenght	240 cm (94,48 in)
Weight	0,05 kg (0,11 lb)
Sterilization	NO
Sterilisation cycles	NO
Compatibility	00900.02, 00900.04

50 D	x	80 D	x	80	x	120	x	160	x	200	x	TOUCH 200	x	FLASH 120	x	FLASH 160	x	FLASH 200	x	300 HP	x	400 HP	x	TOUCH HP	x	TOUCH HP (Argon)	x	SURTRON EVAC	x	SURTRON ABC	x
------	---	------	---	----	---	-----	---	-----	---	-----	---	-----------	---	-----------	---	-----------	---	-----------	---	--------	---	--------	---	----------	---	------------------	---	--------------	---	-------------	---



pcs: 12

00900.FU01/10*



Description	ULPA Filter for SURTRON EVAC - red
Filter Diameter	22 mm 7/8" (0,867 in)
Particles diameter	0,12 micron
Filter Type	ULPA with active carbon
Filter efficiency	99,999X %
Filter duration	00900.FST/10: 5 hours
Filter duration	00900.FU01/10: 20 hours

50 D	x	80 D	x	80	x	120	x	160	x	200	x	TOUCH 200	x	FLASH 120	x	FLASH 160	x	FLASH 200	x	300 HP	x	400 HP	x	TOUCH HP	x	TOUCH HP (Argon)	x	SURTRON EVAC	x	SURTRON ABC	x
------	---	------	---	----	---	-----	---	-----	---	-----	---	-----------	---	-----------	---	-----------	---	-----------	---	--------	---	--------	---	----------	---	------------------	---	--------------	---	-------------	---




pcs: 3



TROLLEY

TR003



Description	Trolley 3 shelvees
Size WxHxD	80x46x46 cm (31x18,1x18,1 in)
Weight	18 kg (34,47 lb)
Shelf size	43x36 cm (16,9x14,17 in)
Color	RAL7035 (Light Grey)
Frame Material	Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.
Wheel type	80 mm pirouette, anterior wheels whit block system.

+ 50 D
 + 80 D
 + 80
 + 120
 + 160
 + 200
 + TOUCH 200
 + FLASH 120
 + FLASH 160
 + FLASH 200
 + 300 HP
 + 400 HP
 + TOUCH HP
 + TOUCH HP (Argon)
 + SURTRON EVAC
 + SURTRON ABC




pcs: 1

WxHxD: 13x47x67 cm

Weight: 18 kg

TR003W



Description	Trolley 3 shelvees
Size WxHxD	80x54x46 cm (31x21,2x18,1 in)
Weight	18,5 kg (34,47 lb)
Shelf size	40x40 cm (15,7x15,74 in)
Color	RAL7035 (Light Grey)
Frame Material	Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.
Wheel type	80 mm pirouette, anterior wheels, whit block system.

+ 50 D
 + 80 D
 + 80
 + 120
 + 160
 + 200
 + TOUCH 200
 + FLASH 120
 + FLASH 160
 + FLASH 200
 + 300 HP
 + 400 HP
 + TOUCH HP
 + TOUCH HP (Argon)
 + SURTRON EVAC
 + SURTRON ABC



pcs: 1

WxHxD: 13x47x67 cm

Weight: 18,5 kg

TR004



Description	Trolley 4 shelvees
Size WxHxD	90x46x46 cm (35,43x18,1x18,1 in)
Weight	20 kg (39,6 lb)
Shelf size	40x36 cm (15,7x14,17 in)
Color	RAL7035 (Light Grey)
Frame Material	Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.
Wheel type	80 mm pirouette, anterior wheels whit block system.

+ 50 D
 + 80 D
 + 80
 + 120
 + 160
 + 200
 + TOUCH 200
 + FLASH 120
 + FLASH 160
 + FLASH 200
 + 300 HP
 + 400 HP
 + TOUCH HP
 + TOUCH HP (Argon)
 + SURTRON EVAC
 + SURTRON ABC



pcs: 1

WxHxD: 13x47x75 cm

Weight: 20 kg

Description Trolley 5 shelvees
Size WxHxD 112x46x46 cm (44x18,1x18,1 in)
Weiht 19 kg (41,9 lb)
Shelf size 43x36 cm (16,9x14,17 in)
Color RAL7035 (Light Grey)
Frame Material Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.
Wheel type 80 mm pirouette, anterior wheels whit block system.

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC

pcs: 1
 shelves WxHxD: 13x47x67 cm
 side columns WxHxD: 104x15x20 cm
 Weight: 27,5 kg



TR005



Description Trolley 5 shelvees
Size WxHxD 112x46x54 cm (44x18,1x21,2 in)
Weiht 24,5 kg (54 lb)
shelf size 43x36 cm (16,9x14,17 in)
Color RAL7035 (Light Grey)
Frame Material Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.
Wheel type 80 mm pirouette, anterior wheels whit block system.

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC

pcs: 1
 shelves WxHxD: 13x47x67 cm
 side columns WxHxD: 104x15x20 cm
 Weight: 27,5 kg



TR005W



Description Support for Cylinder
Size WxHxD 18x6x20 cm (7x2,3x7,87 in)
Weiht 0,5 kg (1,1 lb)
Color RAL7016 (Anthracite Gray)
Frame Material Shelves made with steel, covered with a coat of paint.

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP (Argon)
- + SURTRON EVAC
- + SURTRON ABC

pcs: 1
 WxHxD: 67x47x23 cm



TRCYLSUP01



TRDRAWER01



Description Drawer for trolley type TR W
Size WxHxD 38,5x10x28,7 cm (15x3,9x11,3 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP Argon
- + SURTRON EVAC
- + SURTRON ABC



TRDRAWER01W




Description Drawer for trolley type TR W
Size WxHxD 38,5x10x28,7 cm (15x3,9x11,3 in)

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP Argon
- + SURTRON EVAC
- + SURTRON ABC



TRST001



Description SURTRON TOUCH HP Trolley
Size WxHxD 53x108,5x57,5 cm (20x42,7x22,6 in)
Weight 33 Kg (72,75lb)
shelf size 43,5x10x38 mm (17,x3,9x15 in)
Color RAL7016 (Anthracite gray)
RAL6027 (Light Green)
RAL9005 (Intense Black)
Frame Material Trolley with column made of aluminum with adjustable shelves made with steel, covered with a coat of paint.
Wheel type 80 mm pirouette, anterior wheels whit block system.

- + 50 D
- + 80 D
- + 80
- + 120
- + 160
- + 200
- + TOUCH 200
- + FLASH 120
- + FLASH 160
- + FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP Argon
- + SURTRON EVAC
- + SURTRON ABC



Description
Size WxHxD
Weight
shelf size
Color
Frame Material

Basket for ST Trolley
53x108,5x57,5 cm (20x42,7x22,6 in)
3,2 kg (7 lb)
43,5x10x38 mm (17,x3,9x15 in)
RAL9005 (Intense Black)
Basket made with steel,
covered with a coat of paint.

- x 50 D
- x 80 D
- x 80
- x 120
- x 160
- x 200
- x TOUCH 200
- x FLASH 120
- x FLASH 160
- x FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP Argon
- x SURTRON EVAC
- x SURTRON ABC



pcs: 1

TRSTB01



Description
Size WxHxD
Weight
Color
Frame Material

Support section for ST Trolley
43,5x13,5x38 cm (17x5,3x15 in)
5,75Kg (2,26lb)
RAL6027 (Light Green)
Holder made with steel, covered
with a coat of paint.

- x 50 D
- x 80 D
- x 80
- x 120
- x 160
- x 200
- x TOUCH 200
- x FLASH 120
- x FLASH 160
- x FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP Argon
- x SURTRON EVAC
- x SURTRON ABC



pcs: 1

TRSTP01



Description
Size WxHxD
Weight
Color
Frame Material

Basket for ST Trolley
14,5x30x19,3 cm (5,7x11,8x7,6 in)
2,4 Kg (5,29 lb)
RAL6027 (Light Green)
Basket holder applicable to the
Surtron Touch HP Trolley made with
steel, covered with a coat of paint.

- x 50 D
- x 80 D
- x 80
- x 120
- x 160
- x 200
- x TOUCH 200
- x FLASH 120
- x FLASH 160
- x FLASH 200
- + 300 HP
- + 400 HP
- + TOUCH HP
- + TOUCH HP Argon
- x SURTRON EVAC
- x SURTRON ABC



pcs: 1

TRSTCS1





LED SpA

Via Selciatella, 40
04011 - Aprilia (LT) Italy

www.led.it - info@led.it

Tel. +39 06 92870045
Fax. +39 06 92870046

