

All-inclusive Pricing

- ✓ No subscription required
- ✓ 3-year standard warranty
- ✓ Unlimited users
- ✓ Unlimited cloud storage



Superior Image Quality

- High-quality images that are comparable to cart-based systems
- Imaging settings are optimized based on scanning application
- Optimal results from specialized scanners

Reliable Portability

- Wireless freedom: no cables to trip over or break, and less risk of infections during procedures
- Stable connectivity from our proprietary multi-channel Wi-Fi architecture

Easy to Use

- Point-and-Shoot Ultrasound™ technology removes need to adjust knobs and switches. Artificial intelligence adjusts parameters in real-time.
- Compact & waterproof device that feels like a standard ultrasound scanner



| Model | Frequency | Max Depth | Specialties |
|---------------|-----------|-----------|---|
| C3 HD | 2–6 MHz | 32 cm | Primary Care, EMED, EMS, OB/GYN, Urology |
| L7 HD | 4–13 MHz | 7 cm | Anesthesia, Breast, MSK, Pain, Surgery |
| L15 HD | 6–16 MHz | 5 cm | Anesthesia, Breast, MSK, Surgery |
| C7 HD | 3–10 MHz | 20 cm | Paediatrics, Neonatal, Veterinary |
| PA HD | 1–4 MHz | 40 cm | Critical Care, Cardiac, EMED, EMS, OB/GYN |
| EC7 HD | 3–10 MHz | 15 cm | IVF / Reproductive, OB/GYN, Urology |

Purchase includes:

- M-Mode
- Power Doppler
- Color Doppler
- 1 hour remote in-service training
- 1 hour/month of live telemedicine