
The Gynocular™

Mobile Colposcope for detection of Cervical Cancer

The Gynocular™ is a high quality colposcope in a small package! The Gynocular™ has cutting edge optics with a focal length of 300mm and LED technology. It is lightweight (480g) and easily fits into a pocket or a small bag, making The Gynocular™ the world's first truly mobile colposcope.

Features: Replaceable and rechargeable battery, Red-free / Green filter mode, Handheld or used with a monopod or a tripod, Clinically tested in Sweden, Bangladesh, India and Uganda, CE marked and US FDA approved, Patented technology and trademark.

The Gynocular™ can also be used together with The Gynocular™ Triage to Diagnose (T2D). A smart and user friendly software, working on smartphones and desktops. The T2D enables you to record patient data, review, refer patients, consult, teach and confer colleagues. T2D is also FDA approved for using with The Gynocular™.

Technical Specifications

Magnification (x)	5x	8x	12x
Field of view (mm)	40	30	20
Resolution (lines/mm)	25	40	60
Depth of focus (mm)	9	5	2,5

Benefits with The Gynocular™

- Competes with the top Colposcopes on the market at the fraction of the cost
- It is a optical colposcope with the option to be a video colposcope
- Mirroring capabilities to desktop, laptop or screen
- Mobile, easy to use with high quality optics especially designed for cervical examinations



Manufactured by:


www.gynius.se
info@gynocular.com