

NOVAMAG ED-60

Erectile Dysfunction
Magnetic Therapy

AG 2017 NMT - ED60 17-1



FACTORY:

Uzaycağı Caddesi 1253. Sokak No:5
Ostim - ANKARA / TURKEY
Phone: +90 312 395 0355 • Fax: +90 312 394 6198

sales@novamedtek.com
www.novamedtek.com

NOVAMAG ED-60 is a reforming touch to Vasculogenic Erectile Dysfunction treatment.

The NOVAMAG ED-60 is an advanced innovative solution from Novamedtek that is developed and designed specifically to treat the Male Erectile Dysfunction (Impotence). It is a reforming touch and a state-of-art technology that induces the local controlled angiogenesis for improving the penile functions for patients with vasculogenic erectile dysfunction. The ED-60 maintains and improves the narrowing of the blood vessels of the penis, which is occurred due to peripheral vascular disease, diabetes, cardiovascular disease and others, where the blood vessels supplying and flowing through the penis are narrowed or blocked.


The ED-60 is a painless and anesthesia-free magnetic shockwave therapy device that delivers shocks to the desired area linearly. With its exceptional technological flexible arm and rounded therapy head (8 x 3 cm), the ED-60 provides an exceptional coverage to the whole penile area and indeed to the right and left crura, thus to the muscles that are attaching to the pelvis easily and without positioning.

The extracorporeal painless and safe treatment of the impotence by Novamag ED-60 is applied mainly to 4 main areas including the Penile Shaft, right-left Corpus Cavernosum and the Perineum. A typical treatment is 8 sessions with 2 sessions per weeks, in which 3500-4000 magnetic pulses are applied during 40 minutes. Session duration can be arranged by pulse frequency/min. At the end of the sessions, the linear pulses create a noticed vascular improvement in the angiogenesis, and thus leading to neovascularization or easy blood flow.


The Novamag ED-60 is clinically proven cost-effective impotence treatment equipment that is gifted to patients who are untreated by pharmaceutical solution. It is an innovative solution to the Male sexual trauma, Stress, anxiety and depression that prevents erection.

System Features

- Extracorporeal treatment with exceptional linear head that covers the organ.
- Safe, Painless and Free-anesthesia treatment.
- Non- positioning troubles.
- Touch screen control panel
- Unique Anatomic program.
- No side effects especially for patients who are excluded from medicines due to their medical history.
- Flexible arm that reaches to the angiogenesis areas easily.



Flexible Arm with Linear therapy Head



When prescriptive medicines aren't appropriate, NOVAMAG ED-60 is the solution for improving the male sexual life quality.