



MEDICAL ENDOSCOPES




TISSUE MORCELLATOR

RESECTOSCOPES FOR MORCELLATOR

LASER - RESECTOSCOPES



		Telescope 25°	$\varnothing 4 \times 302$	T3315J
	Working Element use with 550 μm fiber		N5800B-4	
	Working Element use with 1000 μm fiber		N5800B-6	
	Inner Sheath		N5815-2	
	Obturator		N5821A	
	Sheath	F26 $\times$ 180	N5810B-1	
	Obturator	F26	N5820A	
	Fiber clip		P8041A	
	ELLIK Evacuator		P1010	
	Adaptor for ELLIK use with N5815 - 2		P1012F	

All dimensions (mm)

Tissue Shaver



- ☐ Features:
 - Speed controlled by foot switch
- ☐ Technical Data:
 - Max. speed: 1000 r/min, adjustable in 10 steps
 - Suction: 0.08 Mpa
 - Power: 100 Watt
 - Power supply: AC 100 V - 240 V 50/60 Hz
- ☐ Size: 400 mm × 120 mm × 360 mm (W × H × D)
- ☐ Weight: 8 kg

Control Unit YSB - III

N5900E 1



- ☐ Features:
 - Waterproof
 - Speed adjustment

Foot Switch

N5921



- ☐ Features:
 - Brushless motor
 - Waterproof
 - Quick coupler
- ☐ Technical Data:
 - Power of motor: 40 Watt

Hand Piece

N5901




Blade $\varnothing$ 4.5 mm, WL: 400 mm

N5912

OPERATING CYSTOSCOPE

OPERATING CYSTOSCOPE



	Telescope 0° instrument channel F13.5	F19.5 × 191	Y6236-1
	Sheath	F26 × 180	N5810B-1
	Obturator	F26	N5820A
	Sealing Cap	$\varnothing 5 \times 10$	P9703D
	Valve		P9702A