



Videoscopes

### Technical Data and Features

- Integrated light source and image processing
- Compatible with HAWK videoscope
- Automatically identifies videoscope, matches optimal image
- Automatically identifies USB drive, displays remaining capacity
- Brightness / sharpness / saturation adjustment, White balance, Freeze, Photo, Video, Zoom, Square / Circular display
- Power supply: DC 12 V 3.6 A
- Size: 250 x 250 x 95
- Weight: 2 kg

Resolution	Light Source Control	Output	Light Source
1920 x 1080	Manual Automatic brightness adjustment	D V I - D H D M I	20 W L E D

## Video Processor

SD-200B



- Integrated light source and image processing
- Automatically identifies videoscope, matches optimal image
- Compact size, lightweight, easy to carry