


DISPOSABLE ACCESS SHEATH

Disposable Access Sheath



Not available Europe and USA

For Video Uretero-Renoscopes, single-use



Inner/Outer	Length	Packaged	Code No.
F10/12	130	10pc	P10040202
	260		P10040212
	360		P10040222
	400		P10040232
	460		P10040242

Inner/Outer	Length	Packaged	Code No.
F10/12	150	10pc	P10040402
	280		P10040412
	380		P10040422
	420		P10040432
	480		P10040442

Inner/Outer	Length	Packaged	Code No.
F12/14	130	10pc	P10040252
	260		P10040262
	360		P10040272
	400		P10040282
	460		P10040292

Inner/Outer	Length	Packaged	Code No.
F12/14	150	10pc	P10040452
	280		P10040462
	380		P10040472
	420		P10040482
	480		P10040492

Inner/Outer	Length	Packaged	Code No.
F14/16	130	10pc	P100402102
	260		P100402112
	360		P100402122
	400		P100402132
	460		P100402142

Inner/Outer	Length	Packaged	Code No.
F14/16	150	10pc	P100404102
	280		P100404112
	380		P100404122
	420		P100404132
	480		P100404142

Collection Bottle	Capacity	Packaged	Code No.
single-use	300ml	10pc	P10020302

- Low bending resistance, curved sheath-suitable for upper and middle calyces
- Can be direct connected to suction device
- Pressure vent to reduce intrarenal pressure
- Clear surgical field of view - high stone clearance vate
- Hydrophilic coating on sheath and dilator tip for ease of insertion