

Video Cystoscope



- Usage:
 - Control platform for videoscopes
- Features:
 - Combined control unit, light source and monitor 21"
 - Stabilized voltage supply
- Functions:
 - USB 3.0
 - White balance
 - Photo / Video
 - Brightness adjustment + / -
 - Light: On / Off
- Technical Data:
 - Power supply: AC 100 V - 240 V 50/60 Hz
 - 21" LCD monitor
 - Format: 16 : 9
 - Resolution: 1920 x 1080
 - Brightness: 250 cd/m²
 - Contrast ratio: 1000 : 1
 - Video output: HDMI, DVI
 - Power LED light source: 30 W
 - Size: 524 mm x 353 mm x 120 mm (W x H x D)
 - Weight: 12.5 kg

Control Center SD-700BI

P6677EA 1



- Usage:
 - Control platform for videoscopes
- Features:
 - Combined control unit, light source
 - Stabilized voltage supply
- Functions:
 - USB 3.0
 - White balance
 - Photo / Video
 - Brightness adjustment + / -
 - Light: On / Off
- Technical Data:
 - Power supply: AC 100 V - 240 V 50/60 Hz
 - Output voltage: 12V DC 3.8A
 - Power: 45W
 - Lamp: 20W (LED)
 - Resolution: 1920 x 1080
 - Video output: HDMI, DVI
 - Size: 250mm x 250 mm x99 mm (W x H x D)
 - Weight: 2 kg

Control Center SD-200B

P6661E

MEDICAL ENDOSCOPES

Video Cystoscope

Video Processor



OR



Video Cystoscope

Code No.	Insertion Part	Distal Part	Working Length	Working Channel	Field of View	Depth of Field	Deflection
DR060380 11GE	Ø 4.5mm	Ø 4.5mm	380mm	Ø 2.0mm	110°	5-50mm	Up 210° Down 210°