



# High-Definition Musculoskeletal Wireless Ultrasound

## Eliminate Guesswork

Clarius HD3 is the ideal wireless ultrasound scanner for MSK imaging. Make a confident diagnosis with high-definition imaging of muscles, joints, ligaments, tendons, and cartilage both at rest and in motion.



# Diagnose Injuries Instantly with High-Definition MSK Imaging

With 8 times the processing power of most handheld ultrasound systems, Clarius HD3 gives you sharp and clear images of muscle and tissue to assess injuries and perform injections with confidence. With AI assistance on your smart device, getting a great image is easy. Start scanning in seconds, then simply pinch to zoom, slide to change gain, and tap to switch modes.

- Confidently identify joint effusions
- Identify and grade sprains and tears, then monitor progress
- Distinguish between calcific tendonitis and rotator cuff tears
- Identify an Achilles tear and monitor treatment efficacy over time
- Foster patient trust and commitment to your treatment plan
- Visualize your needle to avoid nerves and vessels when injecting



**Clarius L7 HD3  
Linear Scanner**

**Frequency: 4-13 MHz  
Max Depth: 11 cm**




**Clarius L15 HD3  
High-Freq. Scanner**

**Frequency: 5-15 MHz  
Max Depth: 7 cm**




**Clarius C3 HD3  
Curvilinear Scanner**


**Frequency: 2-6 MHz  
Max Depth: 40 cm**




ECU




Supraspinatous




Trv Biceps Tendon



Knee Cartilage



Extensor Tendons Wrist



Long BT



I am thankful that Clarius has brought this level of advanced quality imaging to an ultrasound device that is affordable for clinicians like me. Here's to a future with greater, more relevant tissue diagnostic information to the fingertips of more and more POC clinicians."

**Dr. Greg Fritz, PT, DPT, RMSK - Physiotherapist**

# Superior MSK Imaging in a Highly Portable Scanner



## Clarius Cloud

Manage your exams and build reports anywhere, anytime with unlimited Clarius Cloud storage.



## Easy To Use

Friendly AI-powered app for your Apple or Android automatically delivers the best imaging.



## High-Definition

Octal Beamforming technology and 192 piezoelectric elements deliver high performance and resolution.



## Specialized

AI-powered with specialized workflows for MSK.




## Ultra-Portable

Ideal for small spaces, Clarius fits in your pocket. Take it anywhere!






## Wireless

Steadfast connection. Total freedom. Without cords, cleaning and disinfection are simple and fast.




The more tools that a skilled therapist has in his tool belt, the better it can be. I would definitely consider that to be the mark of somebody who's interested in providing quality service. I couldn't possibly give it a higher recommendation."

**Randy Reinhardt - Patient of Dr. Greg Fritz**




# Advanced MSK Package

The Advanced MSK Package offers more flexibility for users who need additional customizations for examinations and procedures for various joint spaces and MSK related anatomy.



## Dedicated Presets

Optimize imaging with presets for shoulder, elbow, hand/wrist, knee, foot/ankle, hip, spine, and more.



## Specialized Measurements

Save time and improve reporting by having pre-labeled measurements at your disposal.



## Split Screen Mode

Add convenience with the ability to compare left and right MSK anatomy.

## Additional Features

- DICOM
- Needle Enhance
- Elastography

# Clarius Membership

The Advanced MSK Package is included with your Clarius Membership.

## Membership Bundle

With your membership plan, you'll receive a discount on a scanner and the following:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- ✓ Unlimited Clarius Cloud Storage
- ✓ Clarius Live
- ✓ Clarius Education
- ✓ Advanced MSK Package

## One-Time Purchase

You may also purchase your scanner outright with the Advanced MSK Package à la carte:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- 💰 Advanced MSK Package

[REQUEST A QUOTE](#)

Clarius is on a mission to make accurate, easy-to-use, affordable ultrasound available to all medical professionals in every specialty. Our ultra-portable wireless scanners deliver best-in-class imaging with an easy-to-use app for iOS and Android that is powered by artificial intelligence and connected to the cloud. More than two million high-definition scans have been performed using Clarius wireless scanners, now available in over 90 countries worldwide.

## Clarius Mobile Health

2985 Virtual Way #130,  
Vancouver, BC, V5M 4X7  
Canada  
hello@clarius.com  
Phone 1-778-800-9975