




High-Definition Wireless Ultrasound for **Pain Management**

Safe, Precise Procedures

Eliminate the guesswork. With clear visualization of your needle, nerves, and MSK anatomy, Clarius HD3 is the ideal wireless ultrasound system for guiding safe, accurate pain medicine injections.



Visualize Your Needle in High-Definition for Safe Procedures

With 8 times the processing power of most handheld ultrasound systems, Clarius HD3 gives you sharp and clear imaging of nerves and target anatomy to perform more accurate injections with confidence. With greater precision, you can use less medication, improve patient safety, and optimize clinical outcomes, improving the overall patient experience.

- Use ultrasound imaging to plan accurate injections
- Visualize your needle, nerves, and vessels for every procedure
- Target the joint space for painless intra-articular injections
- Impress patients with instant pain relief



**Clarius L7 HD3
Linear Scanner**

**Frequency: 4-13 MHz
Max Depth: 11 cm**




**Clarius L15 HD3
High-Freq. Scanner**

**Frequency: 5-15 MHz
Max Depth: 7 cm**




**Clarius C3 HD3
Curvilinear Scanner**


**Frequency: 2-6 MHz
Max Depth: 40 cm**




Achilles Repair




Patellar Tendon Injection




ACJ Injection



Tarsal Tunnel Tibial Nerve



Tendinosis



Plantar Fascia



Wireless ultrasound enables me to focus more on my patient without having wires get in the way of my nerve block. The third generation Clarius is amazing. The image quality is superior. It connects seamlessly. I mean, it's light, it's small, and it's easier to handle when performing a procedure. It's cheaper, more portable and has greater flexibility than the traditional ultrasound machines."

Dr. David Rosenblum, MD - Interventional Pain Management

Image Clarity You Can Trust for Pain Management



High-Definition

Octal Beamforming technology and 192 piezoelectric elements deliver high performance and resolution.



Easy To Use

Friendly AI-powered app for your Apple or Android automatically delivers the best imaging.



Wireless

Steadfast connection. Total freedom. Without cords, cleaning and disinfection are simple and fast.



Specialized

AI-powered with specialized MSK presets for pain management.




Ultra-Portable

Ideal for small spaces, Clarius fits in your pocket. Take it anywhere!






Cloud and DICOM

Save images to DICOM or manage your exams anywhere, anytime with unlimited Clarius Cloud storage.




It really is a fantastic unit. The image quality of the new HD3 is actually pretty amazing. I mean, I've got to say, I am very impressed. The Clarius wireless performance is phenomenal. I would say that the new Clarius HD3 is almost like a traditional cart-based probe without the wire and the cart. And if you put it behind the iPhone, you literally can't see it. The ergonomics and the size are ideal."

Dr. Alan Hirahara, MD, FRCSC



Advanced MSK Package

The Advanced MSK Package offers more flexibility for users who need additional customizations for examinations and procedures for various joint spaces and MSK related anatomy.



Dedicated Presets

Optimize imaging with presets for shoulder, elbow, hand/wrist, knee, foot/ankle, hip, spine, and more.



Needle Enhance

Enhance the visualization of needles when performing guided procedures.



Specialized Measurements

Save time and improve reporting by having pre-labeled measurements at your disposal.

Additional Features

- DICOM
- Elastography
- Split Screen Mode

Clarius Membership

The Advanced MSK Package is included with your Clarius Membership.

Membership Bundle

With your membership plan, you'll receive a discount on a scanner and the following:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- ✓ Unlimited Clarius Cloud Storage
- ✓ Clarius Live
- ✓ Clarius Education
- ✓ Advanced MSK Package

One-Time Purchase

You may also purchase your scanner outright with the Advanced MSK Package à la carte:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- 💰 Advanced MSK Package

[REQUEST A QUOTE](#)

Clarius is on a mission to make accurate, easy-to-use, affordable ultrasound available to all medical professionals in every specialty. Our ultra-portable wireless scanners deliver best-in-class imaging with an easy-to-use app for iOS and Android that is powered by artificial intelligence and connected to the cloud. More than two million high-definition scans have been performed using Clarius wireless scanners, now available in over 90 countries worldwide.

Clarius Mobile Health

2985 Virtual Way #130,
Vancouver, BC, V5M 4X7
Canada
hello@clarius.com
Phone 1-778-800-9975