


EROND SURG ACCESSORIES

One device for all surgical treatments



The universal radio surgical units with special features - for all kind of surgical procedures.

Standard features:






- Extensive range of accessories
- Adjustable coagulation degree
- Adjustable duration of coagulation
- Trigger exit for the adaption of a smoke evacuation system
- Memory function
- Monitored neutral electrode
- Hand pieces, cables, neutral electrode and electrodes autoclavable
- Large assortment of electrodes - most of the colored shafts are bendable
- Case and foil console are clearly arranged, hygienically designed, easy to use, to clean and to disinfect
- Can also be used in connection with endoscopes.
- All kinds of incisions can be performed and smallest to bigger bleedings are stanchd.



Contents

Category 1- Cutting and Coagulation Electordes	3
Category 2 - Electrode sets.....	8
Category 3 - Special Electrodes ENT	11
Category 4 - Special Electrodes Gynecology.....	14
Category 5 - Special Electrodes for Cosmetic, Plastic & Reconstructive Surgery.....	16
Category 6 - Bipolar Electrodes	18
Category 7 - RF-ReFacing® Accessories	19
Category 8 - Bipolar Forceps	22
Category 10 - Smoke Evacuation System	30
Category 11 - Sterile HF Single Use Electrodes in Packages with 5 pcs.....	32

Category 1- Cutting and Coagulation Electordes


Product / Order No.	Description	Package contents
Multi-Tip Fine Wire Electrode No. 41*	EELBG041	2 electrodes 2 spare wires
	With adjustable wire for depth control Yellow shaft, bendable, angled Wire 0.2 mm (0.007 inch) Spare wire	
	* The Multi-Tip Electrode No. 41 is also part of the Standard Electrode Set No. 300 (see Electrode Sets Category 2)	
Needle Electrodes No.6	EELBG006	2 electrodes
	Yellow shaft, bendable, angled Wire 0.2 mm (0.007 inch)	
Needle Electrodes No.2	EELBG002	2 electrodes
	Yellow shaft, bendable, angled Wire 0.2 mm (0.007 inch)	
Ultrafine, surface coated micro incision needle	Ultrafine needle electrode excellent suitable for epilation and very fine cuts Yellow shaft	
		
Needle electrode	EL2090-2	2 electrodes
	Total length 60 mm, Wire 0.1 mm (0.004 inch,)rigid	

Needle electrode	EL2090-5	5 electrodes
-------------------------	----------	--------------




Total length 60 mm,
 Wire 0.1 mm (0.004 inch,)rigid

Needle electrode angled	EL2090A2	2 electrodes
--------------------------------	----------	--------------



Total length 60 mm,
 Wire 0.1 mm (0.004 inch), rigid

Needle electrode angled	EL2090A5	5 electrodes
--------------------------------	----------	--------------




Total length 60 mm,
 Wire 0.1 mm (0.004 inch), rigid

Loop electrodes


Fine wire 0.25 mm (0.010 inch)
 Yellow shafts, bendable
 Shaft length 38 mm
 Electrodes marked with * are part of the Standard Electrode Set No. 300 (see Electrode Sets Category 2)

Loop electrode No. 15	EELBG015	2 electrodes
------------------------------	----------	--------------




Small, oval loop electrode with half-round tip
 Shaft angled, bendable

Round loop electrode No. 44*	ELROUN44	2 electrodes
-------------------------------------	----------	--------------



Diameter 7 mm, bendable

Round loop electrode No. 24	ELROUN24	2 electrodes
------------------------------------	----------	--------------



Diameter 11 mm, bendable

Round loop electrode No. 43*	ELROUN43	2 electrodes
-------------------------------------	----------	--------------



Diameter 13 mm, bendable

Diamond electrodes	Fine wire 0.25 mm (0.010 inch) Yellow shafts, bendable Shaft length 38 mm Electrodes marked with * are part of the Standard Electrode Set No. 300 (see Electrode Sets Category 2)	
---------------------------	--	--

Diamond electrode No. 53	ELDIAM53	2 electrodes
---------------------------------	----------	--------------



Side length 5 mm

Total diamond length 8 mm, bendable

Diamond electrode No. 45	ELDIAM45	2 electrodes
---------------------------------	----------	--------------



Diamond electrode No. 45*

Side length 7 mm

Total diamond length 11 mm , bendable

Triangle electrodes	Fine wire 0.25 mm (0.010 inch) Yellow shafts, bendable Shaft length 38 mm Electrodes marked with * are part of the Standard Electrode Set No. 300 (see Electrode Sets Category 2)	
----------------------------	--	--

Triangle electrode No. 46*	ELTRIA46	2 electrodes
-----------------------------------	----------	--------------








Side length 7 mm, bendable

Triangle electrode No. 26*	ELTRIA26	2 electrodes
-----------------------------------	----------	--------------



Side length 10 mm, bendable

<p>Rigid needle electrodes</p> <p>Needle Electrode Set - see Electrode Sets</p> <p>Category 2</p>	<p>Inflexible, very sharp, not bendable needle electrodes for quick and precise work, e.g. preparation of fatty tissue.</p> <p>Yellow shafts, angled</p> <p>Total length 45 mm</p> <p>Length of Needle 20 mm</p>
<p>Rigid needle electrode No. 37</p>	<p>EELNAD37 2 electrodes</p>
	<p>Diameter of needle 0.5 mm, rigid</p>
<p>Rigid needle electrode No. 38</p>	<p>EELNAD38 2 electrodes</p>
	<p>Diameter of needle 0.7 mm, rigid</p>
<p>Rigid needle electrode No. 39</p>	<p>EELNAD39 2 electrodes</p>
	<p>Diameter of needle 1 mm, rigid</p>
<p>Blade Electrodes</p> <p>Thin blade for exact working</p>	<p>For Cutting, Coagulation and Preparation</p> <p>Suitable for many surgical applications, e.g. for breast surgeries</p> <p>White shaft, rigid</p>
<p>Short blade electrode No. 05</p>	<p>ELBLAT05 1 electrode</p>
	<p>Total length (incl. shaft) 65 mm, Blade 19 mm, rigid</p>
<p>Long blade electrode No. 10</p>	<p>ELBLAT10 1 electrode</p>
	<p>Total length (incl. shaft) 115 mm, Blade 19 mm, rigid</p>

Coagulation electrodes

Blue shafts

Total length 35 mm

Electrodes marked with * are part of the Standard Electrode Set No. 300 (see Category 2)

Ball electrode No. 48*

ELBALL48

2 electrodes



Ball electrode No. 48*

Shaft bendable

Diameter of ball 2 mm

Ball electrode No. 49*

ELBALL49

2 electrodes



Ball electrode No. 49

Shaft bendable

Diameter of ball 3 mm

Ball electrode No. 47*

ELBALL47

2 electrodes



Rigid, thick needle electrode No. 47*

Diameter of needle 1.5 mm

Length of Needle 10 mm

Shaft angled, rigid

Product / Order No.

Description

Package contents

Ball electrode No. 51*

ELBALL51

2 electrodes




Ball electrode No. 51

Total length 100 mm

Diameter of ball 5 mm, rigid

Category 2 - Electrode sets

Product / Order No.	Description	Package contents
Standard Electrode Set No. 300	ELBU0300	Set of 7 electrodes
		
<p>Cut and Coagulation electrodes</p> <p>All shafts bendable</p> <p>1 Multi-Tip No. 41 - wire 0.2 mm (0.008 inch)</p> <p>1 large loop No. 43 - wire 0.25 mm (0.010 inch)</p> <p>1 small loop No. 44 - wire 0.25 mm (0.010 inch)</p> <p>1 diamond No. 45 - wire 0.25 mm (0.010 inch)</p> <p>1 triangle No. 46 - wire 0.25 mm (0.010 inch)</p> <p>1 thick needle No. 47 diameter 1.5 mm</p> <p>1 ball No. 48 - diameter 2.0 mm</p> <p>Individual shapes of the electrode set are available in packs of 2 each.</p> <p>See Category 1.</p>		
Needle Electrode Set	MEDNADEL	Set of 3 electrodes
<p>Especially suitable for the preparation of fatty tissue</p>		

Rigid needle electrodes

All shafts yellow, bendable

Length of needle 20 mm

1 Needle electrode No.37. Description see Category 1

1 Needle electrode No.38. Description see Category 1

1 Needle electrode No.39. Description see Category 1

Individual shapes of the electrode set are available in



packs of 2 each. See Category 1.

Sclerotherapy Electrode Set

ELSKLESA

Set of 3 electrodes

Used for treatment of teleangiectasia, spider veins and spider nevi.

Uncoated fine wire electrodes

All shafts colour coded, rigid


Total length 35 mm



1 Sclerotherapy electrode No. 9 Description see Category 5



1 Sclerotherapy electrode No. 10 Description see Category 5



1 Sclerotherapy electrode No. 11 Description see Category 5

Individual shapes of the electrode set are available in packs of 2 each. See Category 5.

Finest wire Electrode Set

ELEPISE2

Set of 2 electrodes

Used for treatment of teleangiectasia, spider veins, spider nevi, couperose and for epilation

Coated, insulated finest wire needle electrodes

All shafts colour coded, starr, Total length 35 mm











1 Finest wire electrode No.55.Description see Category 5




1 Finest wire electrode No. 56. Description see Category 5

Individual shapes of the electrode set are available in packs of 2 each. See Category 5.

ENT Electrode Set extra long	ELHNOSET Set of 3 electrodes
	All shafts yellow, bendable 1 ENT needle electrode No. 32
	Description Category 3 1 ENT loop electrode No. 35
	Description Category 3 1 ENT ball electrode No. 34
Description Category 3	
Individual shapes of the electrode set are available in packs of 2 each. See Category 3.	
Varicose Vein Electrode Set for the sclerosis veins	ELVARSET Set of 2 electrodes
	With solid, coated steel core and uncoated, slanted steel tip, which protects the soft tissue. Blue shafts, rigid
	1 Varicose vein electrode No. 57 Diameter of needle 0.3 mm
	Total length of needle 14 mm. Coating 10 mm Total length of electrode 78 mm
	1 Varicose vein electrode No. 58 Diameter of needle 0.35 mm
	Total length of needle 18 mm. Coating 14 mm Total length of electrode 82 mm

Category 3 - Special Electrodes ENT

ENT electrode for Paracentesis	ELHNOP97	1 electrode
--	----------	-------------



Paracentesis electrode No. 97
 for the safe and gentle incision of the eardrum
 Shaft white, rigid
 Needle 50mm, 49 mm coated
 Tip conical

ENT electrodes for nasal radioconchotherapy with partly isolated tips	ELHNOK52	1 electrode
---	----------	-------------



All shafts blue
 Pointed needle
 Total length of needle 26 mm (10 mm insulated), rigid
 Bayonett-shaped ENT electrode No.52 Total length 155 mm

ENT electrodes for nasal radioconchotherapy with partly isolated tips	ELHNOK40	1 electrode
---	----------	-------------



Bent concha electrode No. 40
 according to W. Koenig, MD
 Total length 120 mm. Length from end to bent 45 mm, rigid

Special electrode No. 50 ELHUPP50 1 electrode
 for uvulopalatoplasty



Teflon-coated, for the safe and gentle removal of the uvula and soft palate.
 For the treatment of chronic snorer.
 Yellow shaft, rigid
 Total length 100 mm

Fine wire electrode No. 42 ELCONC42 1 electrode
 for the pharyngeal area



With adjustable wire for depth control
 blue shaft, bendable. Total length 80 mm. Wire 0.2 mm (0.008 inch). Spare wire

ENT electrodes extra long No.32 ELHNON32 2 electrode



Fine wire
 Yellow shaft, bendable
 Total length 100 mm

ENT loop electrode No. 35 ELHNOS35 2 electrode



ENT needle electrode No. 35
 Fine wire
 Yellow shaft, bendable
 Total length 128 mm, Diameter of loop 9 mm

ENT electrode set see category 2 ELHNOK34 1 electrode

ENT ball electrode No. 34



Blue shaft, bendable
 Total length 100mm
 Diameter of ball 2.5 mm

Needle electrode No. 36 ELHNOC36 1 electrode

for coagulation of lingual tonsils and for submucosal UPP



Coated, bent needle electrode
 suitable also for the preparation of fatty tissue
 Blue shaft, rigid
 Sharp needle
 Total length of needle 18 mm (12 mm insulated)
 Total length 105 mm


Category 4 - Special Electrodes Gynecology










BIO-CONE sterile

Rotational Cone Biopsy Electrode

NEW!!! Now available with 4 stabilizer prongs! A very easy method available to perform cone biopsy in an office or clinic under local or general anesthesia with less or no bleeding.

Stabilizing prongs in combination with rotational hub make sure an all over satisfaction in comfort and safety. 9 different sizes are available from small/flat to large/deep



	small / flat	1 electrode	ELBC1812
	small / medium	1 electrode	ELBC1815
	small / deep	1 electrode	ELBC1818
	medium / flat	1 electrode	ELBC2412
	medium / medium	1 electrode	ELBC2415
	medium / deep	1 electrode	ELBC2418
	large / flat	1 electrode	ELBC3009
	large / medium	1 electrode	ELBC3012
	large / deep	1 electrode	ELBC3015

CAUTION! The BIO-CONE electrodes require a handpiece with a shaft diameter of 2.4 mm (3/32"). Handpieces for the ErondSURG can be found in category 9. Conization can be performed only with high frequency units offering a frequency of at least 400 KHz. The higher the frequency the smoother the probe can be extracted.

Gynecological cutting electrodes

Yellow shaft, rigid

Total length 145 mm



Loop 17x15mm 1 electrode ELGYN001



Loop 10x10mm 1 electrode ELGYN002

Loop 15x15mm 1 electrode ELGYN003

Loop 20x15mm 1 electrode ELGYN004




Category 5 - Special Electrodes for Cosmetic, Plastic & Reconstructive Surgery

Disc electrodes for all tangential ablations, e.g. pigmentary disorders. Thin disc for gentle accurate ablation. Applicable for the removal of plain skin lesions: white shafts, bendable.

Small Disc-Elektrode Nr. 80	ELDISC80	2 electrodes
	Total length 40 mm	Disc 5 x 3 mm
Large Disc-Elektrode Nr. 85	ELDISC85	2 electrodes
	Total length 45 mm	Disc 10 x 5 mm

Coated, insulated finest wire needle electrodes Finest wire electrode set. see Category 2

Used for treatment of teleangiectasien, spider veins, spider nevi, couperose and for epilation. The wire is insulated almost to the tip of the electrode which prevents changes of the epidermis. Total length 35 mm

Coated finest wire electrode No. 55	ELEPIL55	2 electrodes
	Yellow shaft	Finest wire 0.2 mm (0.008 inch), rigid
Coated finest wire electrode No. 56	ELEPIL56	2 electrodes
	Red shaft	Fine wire 0.25 mm (0.010 inch), rigid
Uncoated sclerotherapy electrode No. 09	ELSKLE09	2 electrodes
	Yellow shaft, rigid	

Fine wire needle electrodes for radiosclerotherapy. Used for treatment of teleangiectasien, spider veins, spider nevi, couperose and for epilation. Sclerotherapy electrode set see Category 2. Total length 35 mm

Finest wire 0.1 mm (0.004 inch)


Uncoated sclerotherapy electrode No. 10	ELSKLE10	2 electrodes
--	----------	--------------



Red shaft, rigid

Finest wire 0.2 mm (0.008 inch)


Uncoated sclerotherapy electrode No. 11	ELSKLE11	2 electrodes
--	----------	--------------



White shaft, rigid

Finest wire 0.25 mm (0.010 inch)

Needle electrode No. 63 for preparation of fatty tissue	ELPRAE63	1 electrodes
---	----------	--------------



Coated, bent needle electrode



Blue shaft, bendable

Sharp needle




Total length of needle 18 mm (12 mm insulated)




Total length 105 mm

Category 6 - Bipolar Electrodes

<p>Bipolar needle electrode</p> <p>for volume reduction of the nasal concha</p>	<p>ELBICO03</p>	<p>1 electrode</p>
	<p>Bent at the bottom, insulated needles, ergonomic plastic handle, rigid</p> <p>Total length 188 mm</p> <p>needles 16 mm – isolation 6.5 mm</p> <p>Bent after 97 mm – angle 20 degree</p>	
<p>Bipolar needle electrode</p> <p>for volume reduction of the soft palate and the lingual tonsils</p>	<p>ELBISN04</p>	<p>1 electrode</p>
	<p>Bent at the top, isolated needles, ergonomic plastic handle, rigid</p> <p>Total length 193 mm</p> <p>needles 19.5 mm □ isolation 9.5 mm</p> <p>Bent after 169 mm □ angle 5 degree</p>	

Category 7 - RF-ReFacing® Accessories

RF-ReFacing® Ball electrode large	ELFACE08	1 electrode
	<p>Large ball electrode from special alloy</p> <p>Especially applicable for the treatment of face, throat and décolleté</p> <p>White shaft, rigid</p> <p>Diameter of ball 8 mm</p> <p>Total length 30 mm</p>	
RF-ReFacing® Ball electrode small	ELFACE81	1 electrode
	<p>Small ball electrode from special alloy</p> <p>for single, deep wrinkles</p> <p>Bent shaft, white, rigid</p> <p>Diameter of ball 2 mm</p> <p>Total length 25 mm</p>	
RF-ReFacing® Needle electrode blunt	ELFACE82	1 electrode
	<p>Small blunt needle electrode from special alloy for single deep wrinkles</p> <p>Bent shaft, white, rigid</p> <p>Diameter of needle 1.5 mm</p> <p>Length of needle 6 mm</p> <p>Total length 30 mm</p>	

RF-ReFacing® Cone electrode	ELFACE83	1 electrode
 <p>A cone-shaped electrode with a grey shaft and a grey conical tip. The tip is labeled with the number #83.</p>	<p>Cone electrode from special alloy, especially applicable for the treatment of face, throat, décolleté and other larger areas of the body</p> <p>Yellow shaft, rigid. Diameter of cone 15mm. Total length 20 mm</p>	
RF-ReFacing® Plate electrode	ELFACE84	1 electrode
 <p>A plate-shaped electrode with a grey shaft and a black circular plate tip. The tip is labeled with the number #84.</p>	<p>Plate electrode especially applicable for the treatment of larger areas of the body</p> <p>Yellow shaft, rigid</p> <p>Diameter of plate 30 mm</p>	
RF-ReFacing® hand piece	RSWEIS01	1 hand piece
 <p>A long, white, cylindrical hand piece with a tapered end.</p>	<p>Hand piece with 1.6 mm tongs</p> <p>Suitable for the RF-ReFacing® electrodes</p> <p>White colour</p> <p>Length of cable 2.5 m</p> <p>Activation via foot pedal</p>	

Self-adhesive disposable neutral electrode	RS811131	50 neutral electrodes
---	----------	-----------------------



For adults and children (5 kg and more)
 Contact area 105 cm²
 not devided, without cable

Cable for self-adhesive disposable neutral electrodes	RS294030	1 cable
--	----------	---------



Length 4.5 m
 With plug to the unit and double-pinned clamp connection
 Important!
 Please ensure that only Meyer-Haake's disposable neutral electrodes and cables are used. Using other products may lead to problems for which we will not be held liable.

Foot pedal for ErondSURGnot autoclavable	RSFUSS01	1 foot pedal
---	----------	--------------



Length of cable 2.5 m

Holder for hand pieces and electrodes (without accessories)	RFHAMAG1	1 holder
--	----------	----------



Autoclavable holder for hand pieces and electrodes with magnets to affix it to the unit or other metal surfaces

Category 8 - Bipolar Forceps

Bipolar forceps with flat euro plug	RS600520	1 forceps
Standard – short, straight	RS600519	



Short, straight

Total length of forceps 110 mm

Width of tip 0.5 mm

This pair of bipolar forceps is part of the ErondSURG(with bipolar function) basic equipment

Short, straight

Total length of forceps 110 mm

Needle Tip

Bipolar forceps with flat euro plug	RS600521	1 forceps
Short, straight, angled tip	RS600523	



Short, straight, angled tip

Total length of forceps 110 mm

Angled tip

Width of tip 0.5 mm

Short, straight, angled tip

Total length of forceps 110 mm

Angled tip

Needle Tip

Bipolar forceps with flat euro plug Long, straight	1 forceps	Order No.
---	-----------	-----------



Total length of forceps 160 mm. Width of tip 1 mm	RS600539
Total length of forceps 160 mm. Width of tip 2 mm	RS600527
Total length of forceps 160 mm. Needle tip	RS600511
Total length of forceps 195 mm. Width of tip 1 mm	RS600540
Total length of forceps 195 mm. Width of tip 2 mm	RS600529
Total length of forceps 195 mm. Width of tip 2 mm. (Tip with diamond grooves)	RS600557
Total length of forceps 220 mm. Width of tip 1 mm	RS600541
Total length of forceps 220 mm. Width of tip 2 mm	RS600530

Bipolar forceps with flat euro plug Long, straight, angled tip	1 forceps	Order No.
---	-----------	-----------



Total length of forceps 160 mm. Width of tip 0,5 mm	RS600513
Total length of forceps 160 mm Width of tip 1 mm	RS600514
Total length of forceps 160 mm Needle tip	RS600536
Total length of forceps 195 mm Width of tip 1 mm	RS600501
Total length of forceps 195 mm Width of tip 2 mm	RS600507
Total length of forceps 195 mm Width of tip 2 mm (Tip with diamond grooves)	RS600558
Total length of forceps 220 mm Width of tip 1 mm	RS600524

Bipolar forceps with flat euro plug	1 forceps	Order No.
Long, bayonett-shaped		




Total length of forceps 165 mm Width of tip 1 mm	RS600518
Total length of forceps 195 mm Width of tip 1 mm	RS600502
Total length of forceps 195 mm Width of tip 2 mm	RS600534
Total length of forceps 220 mm Width of tip 0.5 mm	RS600531

Bipolar forceps with flat euro plug	1 forceps	Order No.
Long, bayonett-shaped, angled tip		







Total length of forceps 195 mm Width of tip 1 mm	RS600526
Total length of forceps 220 mm Width of tip 1 mm	RS600563

Non-stick bipolar forceps	1 forceps	Order No.
 <p>No sticking of biological tissue</p>	Total length of forceps 110 mm	
	Length of tip 6 mm	RS607020
	Width of tip 0.5 mm	
	Total length of forceps 160 mm	
	Length of tip 8 mm	RS607039
	Width of tip 1 mm	
	Total length of forceps 160 mm	
	Length of tip 8 mm	RS607027
	Width of tip 2 mm	
	Total length of forceps 195 mm	
	Length of tip 8 mm	RS607040
	Width of tip 1 mm	
	Total length of forceps 195 mm	
	Length of tip 8 mm	RS607029
	Width of tip 2 mm	
	Total length of forceps 220 mm	
	Length of tip 8 mm	RS607030
	Width of tip 2 mm	

Non-stick bipolar forceps / angled	1 forceps	Order No.
	Total length of forceps 160 mm	
No sticking of biological tissue	Length of tip 6 mm	RS607014
	Width of tip 1 mm	
	Total length of forceps 195 mm	
	Length of tip 8 mm	RS607001
	Width of tip 1 mm	
	Total length of forceps 195 mm	RS607007
	Length of tip 8 mm	
	Width of tip 2 mm	
	Total length of forceps 195 mm	
	Length of tip 6 mm	RS607002
	Width of tip 1 mm	
	Total length of forceps 220 mm	
	Length of tip 6 mm	
	Width of tip 2 mm	
	Total length of forceps 195 mm	RS607063
	Length of tip 8 mm	
	Width of tip 2 mm	
	No sticking of biological tissue	

Additional Accessories

Neutral electrode	1 piece	RS242003
 	Neutral electrode for ErondSURG Made of conductive silicone rubber 150 x 250 mm 0.5 m connection cable Neutral electrode cable for ErondSURG with flat plug for connection to the neutral electrode Length 3 m	
Punched rubber strip	1 piece	RS194050 RS194100
	With botton for the connection of the neutral electrode Length 0.5 m Length 1 m	
Connection button	1 piece	RS194000
Connection button for rubber strip		
Bipolar cable	1 piece	RSBIPKAB
	Length 2.5 m	
Bipolar cable	1 piece	RSBIPKUS
	for US bipolar forceps Length 2.5 m	

Foot pedal for ErondSURG	1 piece	RSFUSS01
not autoclavable		RSFUSS02



Length of cable 2.5 m

Length of cable 4.0 m

Self-adhesive disposable neutral electrode	50 neutral electrodes	RS811131
not autoclavable		



For adults and children (5 kg and more)

Contact area 105 cm²

not divided, without cable

Cable for self-adhesive disposable	1 piece	RS294030
neutral electrodes		




Length 3.0 m


With plug to the unit and double-pinned clamp connection to the disposable neutral electrode

Important! Please ensure that only Meyer-Haake's disposable neutral electrodes and cables are used. Using other products may lead to problems for which we will not be held liable.

Holder for Hand pieces and	1 piece	RSHAMAG1
bipolar cable (without unit and accessories)		



Autoclavable holder for hand pieces and bipolar cable of the ErondSURG with magnets to affix it to the unit or other metal surfaces

Electrode rack (without electrodes)	ELST0001	1 rack
	<p>Stainless steel jar with cover and removable rack for 12 electrodes with shaft diameter 1.6 mm (1/16").</p> <p>Diameter 80 mm</p> <p>Height 85 mm</p> <p>Jar, cover and rack sterilizable</p>	
Hand pieces for ErondSURG	1 piece	RSGELB01
	<p>Hand piece for cutting for ErondSURG with yellow button</p> <p>Autoclavable - Length of cable 2.5 m</p>	
Hand pieces for ErondSURG	1 piece	RSBLAU01
	<p>Hand piece for coagulation for ErondSURG with blue button</p> <p>Autoclavable - Length of cable 2.5 m</p>	
Special hand pieces for ErondSURG and other radio surgical units	1 piece	RS219030
	<p>Hand piece for ErondSURG for cutting or coagulation - Without finger-switch</p> <p>Must be activated via foot pedal for BIO-CONE (see Category 4) and other electrodes with shaft diameter of 2.4 mm (3/32")</p>	
Length of cable 3 m - RS219030		
Length of cable 4.5 m - RS219045		

Special hand pieces for ErondSURG and other radio surgical units	1 piece	RSH92097
---	---------	----------

Hand piece for cutting or coagulation Without finger-switch - with screw cap Must be activated via foot pedal - for electrodes with shaft diameter of 1.6 mm, Length of cable 3m

Category 10 - Smoke Evacuation System

Smoke Evacuation System smokeSTAR	1 piece	MH100001
--	---------	----------



Smoke evacuation system with pneumatic foot pedal
 230 Volt, 2.5 m power supply connection cable
 Dimensions: 355 x 415 x 205 mm
 Weight (without Filter): 8.4 kg
 Weight (with Filter): 9.1 kg

Trigger Cable	1 piece	MH100006
----------------------	---------	----------

for the adaptation to the radioSURG® automatically activating the smoke evacuation thus granting a long life of the filter
 Length: 2.5 m

4-stage high performance filter	1 piece	MH100002
--	---------	----------








with 3 inputs
 Operational time in normal mode: 32 hours
 Operational time in turbo mode: 16 hours

Vacuum Tube	24 piece	MH100003
--------------------	----------	----------



Tube with suction protection, Length 3.05 m, Ø 22 mm - sterile

Pen Holder	10 piece	MH100004
With adapter for re-usable tube, \varnothing 9.5 mm x length connected to hand pieces on tube, \varnothing 22 mm, 1.83 sterile		
Reducer	12 piece	MH100005
	\varnothing 22 mm to \varnothing 9,5 mm	
Tube	12 piece	MH100007
	Re-usable, \varnothing 22 mm, 3 m	
Speculum Tube Set	5 piece	MH100008
	Tube \varnothing 6.5 mm, 0.61 m to tube \varnothing 22 mm, 1.83 m, sterile	
Extension Tube	24 piece	MH100009
	\varnothing 22 mm 1.83 m	
Remote Control	1 piece	MH100010
	Automatic remote control with 3.05 m connecting cable (if it is not possible to connect a trigger cable)	

Armstand for optimal positioning of the tubes	1 piece	MH100011
Laparoscopy-Tube Set	12 pieces	MH100012



Tube Ø 6.5 mm, 3.05 m, Luer-Connector, valve, filter connector Ø 22 mm, sterile

Category 11 - Sterile HF Single Use Electrodes in Packages with 5 pcs.

Cutting electrode	5 pieces	SELNA02G
--------------------------	----------	----------

Steril with bent shaft and straight wire



Loop electrode 7 mm	5 pieces	SELRO44G
----------------------------	----------	----------

Steril- with straight shaft, loop diameter 7 mm



Loop electrode Set 13 mm	5 pieces	SELRO43G
---------------------------------	----------	----------





Steril- with straight shaft, loop diameter 13 mm



Diamond electrode	5 pieces	SELDI45G
--------------------------	----------	----------

Steril- With straight shaft, side length 7 mm



Triangle electrode	5 pieces	SELTR46G
	Steril - with straight shaft, Side length 7 mm	
Blade electrode	5 piece	SELBL05G
	Steril, with straight shaft, Blade 19 mm	
Ball electrode 2 mm	SELBA48G	5 piece
	Steril with straight shaft, ball diameter 2 mm	
Ball electrode 3 mm	SELBA49G	5 piece
	Steril - with straight shaft, ball diameter 3 mm	