

## MEDICAL FURNITURE



We develop and produce specialized medical furniture for private medical practices. Moreover, we offer you solutions in organizing your private medical practice in order to transform it into a pleasant working place for you and your patients.

Our goal is to create furniture that fits to the chosen environment ergonomically as well as esthetically.




# Negatoscope

Area of use: private medical practices, hospitals, classified as a radiology accessory;

The negatoscope is used to view medical images produced by a chain of high-tech equipment.

## Standard features model NEG-0:

- Dimensions (LxWxH): 90x49x6cm;
- Weight: 4kg;
- Homogeneously illuminated surface;
- Brightness (lux): 0-4000;
- Power supply: 220VAC 50Hz;
- Aluminum frame
- Inox clip
- With remote control
- Mountable on wall



# Germicid lamp

- 3x18W: LG 3x18-M1
- 2x55W: LG 2x55-M1
- 3x55W: LG 165-M1

Germicide lamps are used to air and surface disinfection through UV radiation in spaces with high level of microorganisms eg: hospitals, medical and microbiological laboratories, surgical rooms, private medical practices, waiting rooms, pharmaceutical production halls.

The lamps come with an orientation device that allows the vertical rotation of the lamp with approx. 180° choosing the UV radiation direction.



## Technical description:

<b>UV Tubes</b>	2 or 3 OSRAM UV-C tubes (254nm wavelength) Life last of tubes up to 15% radiation decrease: min. 8000 hours
<b>Input</b>	230V/50Hz-60Hz, with grounding power outlet
<b>Cable length</b>	3m
<b>Other characteristic</b>	Ozone free
<b>Timer</b>	2h
<b>Electric protection class</b>	I, IP20
<b>Warranty</b>	12 months

## Stand and hanging options:

1. Mobile stand with wheels TR-0
2. Fix desktop stand TM-0
3. Wall hanger support TP-0





# Medical examination table

Area of use: private medical practices, emergency rooms;

Medical examination tables are used to comfortably support patients during medical examinations or treatments.

## Standard features model PC-190 TPI-1

- 6cm thick mattress upholstered with eco leather;
- Square inox frame with back rest, adjustable by crank mechanism (PC-190 TPI-1)
- Supporting weight: 120kg
- Dimensions (LxWxH): 190x65x75cm
- Weight: 30kg;



## Optional accessories:

1. Mountable infusion stand on any side of the bed; SS-P-0
2. Paper towel support; SP-P-0
3. Glove support; SM-P-0
4. Leg support: SS-SP-0
5. For larger orders, other leather colors are also available

## Standard colors:



B365\_1044\_AS



B365\_6263\_AS



B365\_9001\_AS

# Medical cart

Area of use: private medical practices, hospitals;

Mobile medical carts keep medical supplies organized, increasing efficiency.

## Standard features model CA-AR-02

- Supporting aluminum column painted in black, on which all the components are mounted;
- Box with silent rail operated and push opening system drawer for accessories and instruments;
- Transparent Plexiglas shelves: 40x49cm sized, height adjustable;
- Screen holder suitable for 15-24" screens;
- 4 power outlets with grounding for medical devices,
- 4m long cable with ON/OFF switch;
- 4 wheels (2 lockable);
- Dimensions (LxWxH): 50x50x137cm;
- Weight: 20kg



## Optional accessories:

- Endoscope and video camera holder with mounting possibility on either side;

## Available colors:

snow	WH10
plum	4H01
pine	8CD1
terra	8CD1



SE-F-0



SE-F-1



SE-R-0

Flexible endoscope holder

Rigid endoscope holder

## Examination bed with drawers

Area of use: private medical practices, emergency rooms;

Medical examination tables are used to comfortably support patients during medical examination or treatments, while the built in drawers and cabinets offer you the necessary space and comfort to store the needed appliances.




### Standard features model PC-190 PAL-0

- 6cm thick mattress upholstered with eco leather;
- Square inox frame with back rest adjustable by crank mechanism
- Supporting weight: 130kg
- Cabinets with 20cm height shelves and jolly basket
- Drawers with silent rail
- Dimensions (LxWxH): 190x65x75cm
- Weight: 100kg;

### Optional accessories:

1. Mountable infusion stand on any side of the bed; SS-P-0
2. Paper towel support; SP-P-0
3. For larger orders, other leather colors are also available
4. Drawers/shelves/baskets at your own choice

### Standard colors:



B365\_1044\_AS



B365\_6263\_AS



B365\_9001\_AS


## Medical cart for radiocauter

Area of use: hospitals, surgical operation rooms;

Mobile medical cart keeps radiocauter and all its accessories organized increasing efficiency

### Standard features model CA-R-0

- Supported aluminum column painted in black
- Box with silent rail and push-open system
- Transparent plexiglass shelf for radiocauter 39x42cm
- Inox shelf for accessories
- Power outlet with grounding for the medical device
- Foot pedal storage place
- 4 wheels (2 lockable)
- Dimensions (LxWxH): 50x50x110cm
- Weight: 14kg



### Optional accessories:

- Endoscope holder with mounting possibility on either side;

### Available colors:

snow	WH10
plum	4H01
pine	9C01
terra	9C01

## Massage bead

Area of use: massage and beauty salons

### Standard features model PC-190 TPI-1:

- 6cm thick mattress upholstered with white eco leather;
- Square inox frame
- Back rest with a face-hole
- Adjustable by crank mechanism
- Supporting weight: 120kg
- Dimensions (LxWxH): 190x65x75cm
- Weight: 30kg;



### Optional accessories:

1. Paper towel support; SP-P-0
2. Glove support; SM-P-0
3. For larger orders, other leather colors are also available

### Standard colors:



B365\_1044\_AS



B365\_6263\_AS



B365\_9001\_AS

## Hospital bassinet

Area of use: hospitals

The hospital bassinet is the first bed of the newborn baby. Newborn medical bassinets are ideal for treating newborn infants.

### Standard features model PNN-0:

- Crib made of non-allergenic nontoxic transparent Plexiglas;
- Trolley made of inox pipe;
- Height adjustable telescopic stand with steel base;
- 4 wheels (2 lockable);
- 5cm thick sponge mattress upholstered with eco leather;
- Dimensions (LxWxH): 85x50x105cm;
- Weight: 3.8kg;



## Changing table

Area of use: private medical practices, hospitals;

A changing table makes diaper changing less of a chore. Drawers and shelves allow you to keep all the necessary items within arms' reach so you're never tempted to leave your baby unattended.

### Standard features model MInf-0

- Made out of 19cm thick plywood;
- Silent rail drawers and high adjustable shelf in a closed cabinet;
- 4 wheels (2 lockable);
- Dimensions (LxWxH): 75x55x80cm;





## Infusion stand

Area of use: private practices, hospitals;

Infusion stands are designed for hanging up infusion liquid jars and infusion pumps. The height is adjustable by screw. All infusion stand models are made of inox bars.

### Standard features:

Supporting weight: 3kg/hook

Base:

- inox 3 wheels (1lockable) SP-200 IP3-1
- silver painted polypropylene 5 wheels; SP-200 SC1-1

Size (HxD): 175cm x 65cm

Tube size:  $\varnothing$ 25x2mm,  $\varnothing$ 20x2mm

Adjustable telescopic stand (H): 80cm-180cm



## Hospital bowl stand

Area of use: private practices, hospitals.  
Made of inox bars.

Base:

- inox 3 wheels (1lockable) B-24-IP3-0
- silver painted polypropylene 5 wheels; B-24- SC1-0

### Standard features model B-24-IP3-0:

Size (DxH): 70x80-110cm

Weight: 3.5kg

Tube size:  $\varnothing$ 25x2mm,  $\varnothing$ 20x2mm

Bowl size:  $\varnothing$ 24cm



## Medical screen

Area of use: private practices, hospitals, beauty salons.

The screen is foldable, equipped with two wheels that allow easy movement to the desired location. Made of tubular inox frame and nonwoven polypropylene material.

### Standard features model P-170-I2-0:

- Open size (LxH): 190x170cm;
- Closed size (LxH): 60x170cm;
- Inox frame pipe size:  $\varnothing$ 20mm,  $\varnothing$ 15mm



SC CATTUS SRL

540236, Targu Mures, Gh. Doja street nr. 177  
Prodcomplex, C1 building  
+4 0724 044 076 | office@erond.ro | www.erond.ro