

EXAMINATION LED LAMPS

MIMSAL



MODEL (code)	MS LED PLUS (M187001)	LUXIFLEX LED PLUS (I0104NL)	MS LED (M183500)	LUXIFLEX LED (I0104AOLED)	LS LED (LS010316LED)
Light source	LED	LED	LED	LED	LED
Watt	12W	10W	8W	8W	7,5W
Lux at 50 cm	45.000 lux	35.000 lux	25.000 lux	15.000 lux	2.000 lux
Colour temperature	4.500 °k	4.500 °K	4.500 °k	4.000 °K	4.000 °K
Colour rendering	97%	97%	97%	80%	80%
Lamp lifetime	50.000 h	50.000 h	50.000 h	25.000 h	15.000 h
Arm length	90 cm.	72 cm.	90 cm.	72 cm.	90 cm.