

KATRACORE® AUTOMATIC BIOPSY NEEDLE

- ↳ Easy and one hand operated with two stage arming
- ↳ Echogenic and atraumatic needle tip design
- ↳ 22mm advancement with 20mm notch
- ↳ Ergonomic and light weight design
- ↳ Depth marked canula
- ↳ Coaxial needles are optionally available



Ref. No	Diameter (G)	Length (cm)	Coaxial Needle Length (cm)	With Coaxial Ref. No
GKR1410	14	10	5,5	GKR1410-X
GKR1415	14	15	10,5	GKR1415-X
GKR1420	14	20	15,5	GKR1420-X
GKR1610	16	10	5,5	GKR1610-X
GKR1615	16	15	10,5	GKR1615-X
GKR1620	16	20	15,5	GKR1620-X
GKR1810	18	10	5,5	GKR1810-X
GKR1815	18	15	10,5	GKR1815-X
GKR1820	18	20	15,5	GKR1820-X
GKR1825	18	25		
GKR1830	18	30		
GKR2010	20	10	5,5	GKR2010-X
GKR2015	20	15	10,5	GKR2015-X
GKR2020	20	20	15,5	GKR2020-X

