



CARING
FOR WOMEN FOR
OVER 18 YEARS



- SONY Camera
- 40x Magnification
- 3 Green Filters
- Plug and Play
- Preinstalled Workstations Available



Let Us Be Your


Partner In Advancing Medical Care.



105 Remington Blvd, Suite C, Ronkonkoma, NY 11779, U.S.A
www.lutechmedical.com

Ph: 631-676-7432 | Fax: 631-619-0723
info@lutechmedical.com

COLPOSCOPE



High color rendering and low color temperature produce an increase in the contrast of the cervical tissue imaging result. Improved imaging results allow for more accurate diagnosis and treatment.

LT-300 Digital Video Colposcope

The LT-300 is a class leading Digital Video Colposcope System that is fitted with a high definition Sony CCD imaging module, and an astounding 550 lines of horizontal definition. Boasting a super bright LED cold light source, the LT-300 is able to display true life colors of the tissues. Equipped with a three grade green electronic filter, the LT-300 features a acetic acid timing and magnification display, and is able to render 1-36x optical magnification and 1-40 x magnification. There are a wide range of focus positions to choose from to include micro, manual and auto positioning, as well as an optional intuitive management software.

Integrate image and video into your current Patient Management Technology.

- SONY High Definition CCD Imaging Module
- 40X Digital Magnification
- Three Grade Electronic Green Filters
- Auto focus and Manual focus modes
- 550 TV lines. Available in PAL and NTSC
- Magnification and Acetic acid timing display
- Super Bright LED light source
- Ergonomic Design
- Plug and Play connection to your computer
- Easily move images to current software
- Available with optional pre installed workstation computer
- Optional Software Package Available

