

# MS 1400



## Electronic Gynecology Table



### MESPA HEALTHCARE PRODUCTS INC.

2. organize sanayi bölgesi 83216 nolu cadde no:20 ŞehitkamilGaziantep / TURKIYE

Tel: +90 444 5 821

Fax: +90 342 337 53 13

e-mail: [export@mespa.com.tr](mailto:export@mespa.com.tr)

web: [www.mespa.com.tr](http://www.mespa.com.tr)



# MS 1400



## Electronic Gynecology Table

- 3 motors
- Electronic adjustments of backrest, height
- Trendelenburg functions
- Control Pedal
- Stainless steel main-frame
- Detachable footrest part
- Removable and adjustable Polyurethane leg supports
- Synthetic leather coated mattress
- Stainless steel handles
- Stainless steel waste bin



### TECHNICAL DETAILS

External Dimensions:	75x120cm
Mattress Dimensions:	60x160cm
Height Range:	60-85cm
Trendelenburg Angle:	0-12°
Backrest Angle:	0°-70°
Safe Working Load Capacity:	180kg

### ELECTRICAL DATA

Power In:	max. 3,15A 100V – 240V ~ 50/60 Hz
Electrical Shock Protection:	Class I, Type B
Liquid Ingress Protection:	IPx66

### Color Options:

D101	D103	D106	D107	D108
Claret	Red	Beige	Dark Blue	Black

## MESPA HEALTHCARE PRODUCTS INC.

2. organize sanayi bolgesi 83216 nolu cadde no:20 ŞehitkamilGaziantep / TURKIYE

Tel: +90 444 5 821

Fax: +90 342 337 53 13

e-mail: [export@mespa.com.tr](mailto:export@mespa.com.tr)

web: [www.mespa.com.tr](http://www.mespa.com.tr)

