

MedGyn HSG Catheter

MedGyn's HSG catheter are used for administering contrast media or saline during Hysterosalpingography and Sonohysterography procedures to detect uterine pathology and tubal occlusion. MedGyn's HSG catheter is a disposable balloon catheter with a rigid insertion sheath. A syringe is included for balloon inflation.

Features & Benefits

- Ideal for hysterosalpingography and sonohysterography
- Sterile, disposable, and latex free
- 5Fr. and 7Fr. sizes for easy insertion
- Placement sheath decreases the need for a tenaculum
- Allows for clear view of external cervical OS

022726 HSG Catheter 5Fr (10/bx)

022728 HSG Catheter 7Fr (10/bx)



TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

t: +1 630.627.4105 toll-free: 800.451.9667

f: +1 630.627.0127

e: info@medgyn.com

w: medgyn.com

HSGCATHR-PF2

