


Tumescentor™ System

Tumescent Infiltration for
Endovenous Thermal Ablation



Features

- Operation status display monitor
- Injection speed control
- Dedicated IV pole included
- Tumescentor™ accessories are compatible with:
 - RF Medical RFP -300 pump (hand piece option only)
 - Nouvag pump (Simple tubing & hand piece option)

Specifications

Protection Class	I
Protection Type	BF
Input Power Voltage	AC 100~240V
Input Power Frequency	50/60Hz
Maximum Input Power	100VA
Display	20×2 Character LCD
RPM	30~400RPM
Injection Control	By Foot Switch

Tumescentor™ System



Specifications

Product Code*	Compatible Pump	Hand piece / simple tubing	Injection Needle	Tube Length
TI-1B4PP-2380	RF Medical RFP-300	Hand piece	23G x 80mm	3m
TI-1BNV-3MH	Nouvag series DP-30	Hand piece	N/A	3m
TI-1BNV-4MH				4m
TI-1BNV-3M		Simple tubing		3m
TI-1BNV-4M				4m

* More product codes are available upon request with MOQ.