




Tumescent Infiltration Pump RFP-300

Features

- Operation status display
- Injection speed control
- Dedicated IV Pole included



Specifications

Protection Class	I
Protection Type	BF
Input Power Voltage	AC100V ~ 240V
Input Power Frequency	50/60Hz
Maximum Input Power	60VA
Display	20*2 Character LCD
RPM	30 ~ 400RPM
Injection Control	By Foot Switch

*Specifications are subject to change without notice.



Tumescent Infiltration Hand-Pieces



<Tumescent Infiltration Hand-Piece>

Features

- ⦿ Time saving
- ⦿ Ergonomic designs
- ⦿ Wider selection

Specifications

Type	Product Code	Code Reading	Remark
TI-2B	"TI-◆◆○●-□□△△ [TI-2B○●-□□△△] * Needle Gauge (17~ 30G) * Needle Length (13~190mm)"	◆◆ : Number of Bags (2B: 2Bag) ○ : Pump tube OD size (4 : Φ4.0mm~ Φ4.9mm) (5 : Φ5.0mm~ Φ5.9mm) (6 : Φ6.0mm~ Φ6.9mm) ● : Spike connection type (P : Removable pump tube and spike) (S : Straight connection type) (A : All connection type) □□ : Needle Gauge (17~ 30G) △△ : Needle Length (13~140mm)	
TI-1B4S	"TI-◆◆○●-□□△△ [TI-1B4S-□□△△]"		
TI-1B4P	"TI-◆◆○●-□□△△ [TI-1B4P-□□△△]"		
TI-1B4A	"TI-◆◆○●-□□△△ [TI-1B4A-□□△△]"		
TI-1B4PP-XX	"TI-◆◆○●-XX [TI-1B4PP-XX]"	◆◆ : Number of Bags (1B : 1Bag) ○ : Pump tube OD size (4 : Φ4.0mm~ Φ4.9mm) ● : Spike connection type (PP : Removable pump tube and spike AND Injection Tube ID is Φ1.2mm) XX : Non Injection Needle type □□ : Needle Gauge (21~24G) △△ : Needle Length (60~80mm)	
TI-1B4PP	"TI-◆◆○●-□□△△ [TI-1B4PP-□□△△]"		
TI-1B4NV	"TI-◆◆○○●- [TI-1B4NV-□□△△]"	◆◆ : Number of Bags (1B : 1Bag; 2Bag) ○ : Pump tube OD size (4 : Φ4.0~4.9) ●● : Needle connection type (NV : Compatible type with other company's product) □□ : Injection tube length (3M : 3m, 4M : 4m) △ : have the handle or non-handle (H : have the handle)	1. Non-handle Type 2. Handle Type