



VeinCLEAR™ System

Endovenous Radiofrequency Ablation System

CE
1639



VeinCLEAR™ System Catheter



Heating Element Length



VeinCLEAR™ System offers

- Outpatient procedure
- Minimum complications
- High occlusion rate
- Short recovery time

Features

- 5Fr and 7Fr catheters
- Heating element length from 1cm to 7cm
- Available in Alloy52, Stainless Steel, LED and Guidewire options

REF RCT-003(E) Rev. 11 (2022.03.28)

EC REP

Obelis S.A.
Bd Generali Waris, 53, B-1030
Brussels, Belgium
Tel: +32 (0) 732-5954 Fax: +32 (0) 732-6003
Web: www.obelis.net



RF Medical Co., Ltd.
#502-507, #5T1, #601, World Meridian 254, Beotkkot-ro,
Geumcheon-gu, Seoul, South Korea
Tel: +82 (0) 2108-4200 Fax: +82 (0) 213-7999
E-mail: overseas@rfaco.kr Web: www.rfaco.kr

01 VeinCLEAR™ System

Dedicated for Endovenous Thermal Ablation

+ Generator V-700

- LCD touch screen
- Compatible with all STS and Alloy52 catheters
- Advanced automatic self-test system
- Safe temperature control algorithm
- PC compatibility for observation and data download
- Multi-language support with audio guidance
- Target temperature and ablation cycle time can be optionally adjusted
- Pre-set ablation algorithm available according to vein diameter




Specifications

Product Code	V-700
Protection Class	I
Protection Type	BF
Input Power Voltage	100 - 240V (AC)
Input Power Frequency	50/60 Hz
Maximum Input Power	150 VA
Measuring Temperature	0°C - 200°C
Alarm Sound	More than 65 dB

**M-3004 and V-1000 generators are compatible with Alloy52 catheters only.

03 VeinCLEAR™ Catheters

Dedicated for Endovenous Thermal Ablation



+ Options



Specifications

Product Code	Compatible Introducer	Heating Element	Total Length	Heating Element Length	Compatible Guidewire
VWT 70FS	5Fr	Stainless Steel	70cm	7cm	N/A
VWT 100FS			100cm		
VWT 70FG			70cm		
VWT 100FG			100cm		
VWT 70FGL*			70cm		
VWT 100FGL*			100cm		
VWT 70FGW**	7Fr	Alloy52	70cm	3cm	0.025"
VWT 100FGW**			100cm		
VWT 4510			1cm		
VWT 4515			1.5cm		
VWT 4520			2cm		
VWT 4530			3cm		
VWT 60			7cm		
VWT 7020			2cm		
VWT 7030			70cm		
VWT 70			3cm		
VWT 90	7cm				
VWT 100	100cm	7cm			

* Lead light on tip of catheter ** Embodied guidewire on tip

Compatibility of catheters with generators

	V-700	M-3004	V-1000
STS Catheter	LED	X	X
	General(5Fr & 7Fr) & Guide	O	X
Alloy Catheter	General	O	O

02 Tumescer™ Infiltration Handpieces

Tumescer Infiltration for Endovenous Thermal Ablation

+ Tumescer™

- Tumescer™ enables to decrease the tumescer injection time significantly, and eventually the overall procedure time
- Tumescer™ can be offered with various options



Specifications

Product Code	Type	Injection Needle	Tube Length
TI-IB4PP-2380	Handpiece	23G x 80mm	3m
TI-IBNV-3MH	Handpiece	N/A	3m
TI-IBNV-4MH			4m
TI-IBNV-3M	Simple tubing	N/A	3m
TI-IBNV-4M			4m