

LS

LS is a examination light with articulated arm and metal structure for medical offices.

Easy to use

Adjustable height articulated arm, allows easy adjustment in different positions, provides a long reach and maximum vertical and horizontal movement, and always remains in the correct position.

Switch located on the connection cable.

Available in two versions

LS includes a fluorescent light source.

LS LED includes a LED light source.

Simple and durable

Metallic structure, manufactured with durable materials and with very low maintenance costs.

Mounting Options

Standard version includes AH table clamp.

As accessories, we offer a wide selection of supports such as wall clamps, wall extension arm and trolley stand.




LS

Technical Details

Light source	Fluorescent 15W
	LED 7,5W
Head diameter	21 cm
Socket	E-27 ceramic
Arm length	90 cm
Colour	White

LS versions:

- **LS** (Cod. LS010316)
- **LS LED** (Cod. LS010316LED)



Wall support
Cod. B0101800

AH table clamp
Cod. B010190

Rail Plus bracket
Cod. 82026

Extension wall arm
Cod. B0101S0

Trolley stand
4 kg Cod. 09655