

Video Cystoscope HD

Insertion Part: \varnothing 5.2 mm
 Distal Part: \varnothing 5.2 mm
 Working Length: 380 mm
 Working Channel: 2.0 mm
 Field of View: 100°
 Depth of Field: 3 - 50 mm
 Deflection: Up 210°, Down 120°



Video Uretero - Renoscope HD L - Type

Insertion Part: Fr 8.5
 Distal Part: Fr 7.5
 Working Length: 670 mm
 Working Channel: Fr 3.6
 Field of View: 110°
 Depth of Field: 2-50 mm
 Deflection: Up 270°, Down 270°
 8 h Continuous use



Control Center SD - 300A



HD Camera Head for
Rigid Telescopes



Compact Video Center
SD - 700B

1080P
 7 Modes for rigid endoscopes
 BCR: Image enhancement to reduce
 blood color
 BNB: Image enhancement of blood
 vessels and mucosa structure

1080P
 21" Monitor
 7 Modes for rigid endoscopes
 BCR: Image enhancement to
 reduce blood color

Video Uretero - Renoscope L - Type System

